

Acid Purification System and Cleaning System



Acid Purification System(APS)



Cleaning System(CS)

Features of Acid Purification System(APS) and Cleaning System(CS)

1. Made with Injection molding PFA Teflon.
2. Consist of All PFA and PTFE Supporter.(Metal Free)
3. Heating system with graphite heater.
4. Simplified Installation.
5. Acid Purification System is compatible with Cleaning system.
6. Connectable each other.(refer to below images)



Acid Purification System to Acid Purification System



Acid Purification System to Cleaning System

Features of Acid Purification System

1. High-productivity. (produce over 1.4ml per minute)
2. Acid Purification System for Nitric acid, Hydrochloric acid, Hydrofluoric acid, etc.
3. Easy to clean inside of APS for first use.
4. Compatible with Cleaning system.

Features of Cleaning System

1. Easy to clean Vials, Vessels, ICP, ICP-MS consumables.
2. Clean with Humic acid vapor.
3. Clean with Dilute acid. (by boiling)
4. Recycling acid. (Low maintenance cost)
5. Easy to clean small components.