


인증표준물질

Certified Reference Materials





중금속분석을 위한 시료전처리장치와 인증표준물질등 소모품 전문기업!

(주) 오디랩은 2008년 8월에 설립된 회사로 중금속분석에 사용되는 흑연블럭 산 분해장비인 에코프리 I, II, III 시리즈와 산 세척장치, 고순도 산 제조 장치, 유리분주기, ICP/ICP MS 소모품, 인증표준물질(CRM), 숙련도 평가물질 (PT) 등 을 제조, 수입판매하고 있습니다.

(주) 오디랩에서 제조 판매하는 흑연블럭 산 분해장비는 열선 가열판이나 마이크로웨이브의 단점을 보완한 제품으로 국내를 비롯하여 세계 7개국에 특허를 획득하였고 현재 해외로도 수출 중 에 있습니다.

또한 실험실에서 분석 데이터의 신뢰성확보를 위한 인증표준물질(CRM)과 표준물질(RM), 국제숙련도 물질을 전세계에서 수입하여 판매하고 있습니다. 인증표준물질은 고객이 찾으시는 제품을 탐색하여 드리고 있으며, 가장 근접한 제품으로 추천드리고 있습니다.

특히 유럽환경규제인 RoHS에 대응한 IEC62321시험법에 나오는 인증표준 물질을 국내 시험평가기관이나 국가기관에 공급하고 있으며, 환경부에서 실시하는 정도관리에 대응하여 LGC사에서 제공하는 환경관련 숙련도 물질을 공급하고 있습니다.

**저희 (주) 오디랩은 화학실험실의 동반자로서
분석의 재현성과 정확성, 신뢰성 확보를 위해
언제나 고객의 노력과 함께 하겠습니다**




ODLAB

자동 산분해장비

ADS25



견적문의




경량화 & 벤틸레이션

경량화 & 벤틸레이션 -



산순환 포집분해용기

산순환 포집분해용기 -



오염방지&내구성

오염방지&내구성 -



앱 연동 조작


앱 연동 조작 -

이동 및 보관을 위해 플레이트를 타공디자인으로 경량화를 하였습니다 또한 타공을 통하여 원활하게 열기의 순환이 이루어 집니다.

좌우에 있는 리프트 장치로 산 순환 포집분해 용기 내부의 산을 가열 / 냉각 시켜서 사이펀 현상에 의해 리사이클 시켜 시료를 분해할 수 있도록 디자인 하였습니다.


SUS 재질사용 및 테프론 특수코팅을 하여 표면이 쉽게 오염되지 않도록 제작하였습니다. 또한 상부 가열부와 하부 전자제어부는 서로 격리, 밀봉되어 열 또는 산 증기로 인하여 전자제어장치가 손상되지않도록 되어있습니다.

20 Step 으로 가열 / 냉각 으로 분해조건을 프로그램화 할 수 있으며, 앱을 사용하여 조절가능합니다.



컴팩트한 사이즈

컴팩트한 사이즈 -



오토메틱 리프팅

오토메틱 리프팅 -

메뉴얼 및 프로그래밍 기능으로 반복적인 가열 / 냉각을 할 수 있도록 리프팅 기능이 있습니다.



균일한 온도

균일한 온도 -

흑연 소재를 사용하여 균일한 온도를 제공하고 ($\pm 1^{\circ}\text{C}$ 온도편차를 갖는다) 제어는 0.2°C 로 제어된다.



수동 승강버튼

수동 승강버튼 -

리프트 장치를 수동버튼을 사용하여 상부 랙(Rack)을 상하로 움직여 사용자가 원할 시 용기의 상태를 언제든지 확인 할 수 있습니다.



GLASS EXPANSION
Quality By Design

ICP-OES / ICP-MS

모든 메이커 (애질런트, 씨모, 퍼킨..etc) 전제품



견적문의

PESTICIDE MIXTURES



California Pesticide Class 1 Mixture cis and trans Chlordane

DRE-A50000078AL 1ml
California Pesticide Class 1 Mixture cis and trans Chlordane
100 µg/mL in Acetonitrile(‡)(*)

DRE-A50000079AL 1ml
California Pesticide Class 1 Mixture cis and trans Chlordane
100 µg/mL in Acetonitrile Second Source (‡)(*)

DRE-S50000078AL 5x1ml
California Pesticide Class 1 Mixture cis and trans Chlordane
100 µg/mL in Acetonitrile(‡)(*)

Aldicarb	Carbofuran
Cis-Chlordane (Alpha Isomer)	Coumaphos
Dimethoate	Ethoprophos
Fipronil	Imazalil
Paclbutrazol	Parathion-methyl
Thiacloprid	Trans-Chlordane (Gamma Isomer)
Chlorfenapyr	Chlorpyrifos
Daminozide	Dichlorvos
Etofenprox	Fenoxycarb
Methiocarb	Mevinphos
Propoxur	Spiroxamine

California Pesticide Mixture 1-100

DRE-GA09000471AL 1ml
California Pesticide Mixture 1-100 100 µg/mL in Acetonitrile(‡)(*)

abamectin	acephate
aldicarb	azoxystrobin
bifenthrin	boscalid
chlorantraniliprole	chlorfenapyr
cypermethrin (mix of isomers)	daminozide
acequinocyl	acetamiprid
baythroid (mixture of isomers)	bifenazate
carbaryl	carbofuran
chlorpyrifos	clofentezine
diazinon	dichlorvos

California Pesticides Class 1 Mixture

DRE-GA09000667AL 1ml
California Pesticide Class 1 Mixture
100 µg/mL in Acetonitrile(‡)(*)

DRE-GS09000667AL 5x1ml
California Pesticide Class 1 Mixture
100 µg/mL in Acetonitrile(‡)(*)

DRE-GA09001033AL 1ml
California Pesticides Class 1 Mixture
100 µg/mL in Acetonitrile Second Source(‡)(*)

aldicarb	carbofuran
chlorpyrifos	coumaphos
dimethoate	ethofenprox
fipronil	imazalil
paclbutrazol (mixture of stereo isomers)	phosdrinTM (mevinphos)
thiacloprid	spiroxamine (mix of isomers)
chlordane (mix of isomers)	chlorfenapyr
daminozide	dichlorvos
ethoprophos (prophos)	fenoxycarb
methiocarb	methyl parathion
propoxur	

California Pesticides Class 2A Mixture

DRE-GA09000668AL 1ml
California Pesticide Class 2A Mixture
100 µg/mL in Acetonitrile(‡)

DRE-GS09000668AL 5x1ml
California Pesticide Class 2A Mixture
100 µg/mL in Acetonitrile(‡)

DRE-GA09001034AL 1ml
California Pesticides Class 2A Mixture
100 µg/mL in Acetonitrile Second Source(‡)

abamectin	acephate
azoxystrobin	baythroid (mixt. of 4 isomers)
boscalid	captan
clofentezine	cypermethrin (mix of isomers)
etoxazole	fenhexamid
fludioxonil	fonicamid
acequinocyl	acetamiprid
bifenazate	bifenthrin
carbaryl	chlorantraniliprole
diazinon	dimethomorph (two isomers)
fenpyroximate (as raceimers)	

California Pesticides Class 2B Mixture

DRE-GA09000669AL 1ml
California Pesticide Class 2B Mixture
100 µg/mL in Acetonitrile(‡)(*)

DRE-GS09000669AL 5x1ml
California Pesticide Class 2B Mixture
100 µg/mL in Acetonitrile(‡)(*)

DRE-GA09001035AL 1ml
California Pesticides Class 2B Mixture
100 µg/mL in Acetonitrile Second Source(‡)(*)

dibrom	hexythiazox
malathion	metalaxyl
pentachloronitrobenzene	permethrin (mix of isomers)
prallethrin	propiconazol (mix of isomers)
spinetoram (mix of isomers)	spinosad (mix of spinosyn A & D)
Systhane TM	tebuconazol (Folicur)
imidacloprid	kresoxim methyl
methomyl	oxamyl
phosmet	piperonyl butoxide
pyrethrin (mix of isomers)	pyridaben
spiromesifen	spirotetramat
thiamethoxam	trifloxystrobin

California Supplemental Cannabis Pesticide Mixture 463

DRE-GA09000463AL 1ml
California Supplemental Cannabis Pesticide Mixture 463
100 µg/mL in Acetonitrile(‡)(*)

DRE-GS09000463AL 5x1ml
California Supplemental Cannabis Pesticide Mixture 463 100
µg/mL in Acetonitrile(‡)(*)

captan
dimethomorph
pentachloronitrobenzene
spinetoram (mix of isomers)
coumaphos
fenhexamid
phosdrin TM (mevinphos)

Canada Pesticide Mixture 1

DRE-GA09001041AL 1ml
Canada Pesticide Mixture 1 50 µg/mL in Acetonitrile(±)(*)

abamectin (mix of isomers)	acetamiprid
boscalid	carbofuran
diazinon	dichlorvos
ethoprophos (prophos)	fenpyroximate (raceimers)
malathion	metalaxyl
novaluron	oxamyl
piperonyl butoxide	propoxur
spiromesifen	spirotetramat
thiacloprid	thiamethoxam
aldicarb	bifenazate
chlorantraniliprole	daminozide
dimethoate	dinotefuran
flonicamid	imidacloprid
methiocarb	methomyl
paclobutrazol (stereo isomers)	phosmet
Spinetoram (spinetoram J & L)	spinosad (Mix of Spinosyn A & D)
Systhane TM	tebuconazole
thiophanate methyl	

Canada Pesticide Mixture 2 ver. 2

DRE-A50000072AL 1ml
Canada Pesticide Mixture 2 ver. 2 20-1000 µg/mL in Acetonitrile(±)
DRE-S50000072AL 5x1ml
Canada Pesticide Mixture 2 ver. 2 20-1000 µg/mL in Acetonitrile(±)

Abamectin [100 µg/mL]	Acephate [20 µg/mL]
Chlorpyrifos [40 µg/mL]	Clofentezine [20 µg/mL]
Diazinon [20 µg/mL]	Dichlorvos [100 µg/mL]
Fensulfothion [20 µg/mL]	Fenthion [20 µg/mL]
Metalaxyl [20 µg/mL]	Methomyl [50 µg/mL]
Prallethrin [50 µg/mL]	Propiconazole [100 µg/mL]
Pyridaben [50 µg/mL]	Resmethrin [100 µg/mL]
Tetramethrin [100 µg/mL]	Teflubenzuron [50 µg/mL]
Allethrin [200 µg/mL]	Bifenthrin [1000 µg/mL]
Coumaphos [20 µg/mL]	Cypermethrin [300 µg/mL]
Dimethoate [20 µg/mL]	Ethoprophos [20 µg/mL]
Fipronil [60 µg/mL]	Kresoxim-methyl [20 µg/mL]
Paclobutrazol [20 µg/mL]	Phenothrin [50 µg/mL]
Propoxur [20 µg/mL]	Pyraclostrobin [20 µg/mL]
Spirotetramat [20 µg/mL]	

Canada Pesticide Mixture 2B

DRE-GA09001039AL 1ml
Canada Pesticide Mixture 2B 100 µg/mL in Acetonitrile(±)(*)
DRE-GS09001040AL 5x1ml
Canada Pesticide Mixture 2B 100 µg/mL in Acetonitrile(±)(*)

dibrom
fenthion
methyl parathion
spiroxamine (mix of isomers)
dimethomorph (GC1: 56.8%, GC2: 42.9%)
imazalil
spirodiclofen
tetrachlorvinphos (ISO)

Canada Pesticide Mixture 3 ver. 2

DRE-S50000073AL 5x1ml
Canada Pesticide Mixture 3 ver. 2 20-1000 µg/mL in Acetonitrile(±)(*)

Azadirachtin A [1000 µg/mL]	Chlorfenapyr [50 µg/mL]
Clothianidin [50 µg/mL]	Cyantraniliprole [20 µg/mL]
Daminozide [100 µg/mL]	Dodemorph [50 µg/mL]
Etridiazole [30 µg/mL]	Fludioxonil [20 µg/mL]
Fluopyram [20 µg/mL]	MGK 264 isomer A [50 µg/mL]
Naled [100 µg/mL]	Parathion-methyl [50 µg/mL]
Pyrethrins [50 µg/mL]	

Canada Pesticide Mixture 1 ver. 2

DRE-A50000070AL 1ml
Canada Pesticide Mixture 1 ver. 2 10-1000 µg/mL in Acetonitrile(±)(*)
DRE-S50000070AL 5x1ml
Canada Pesticide Mixture 1 ver. 2 10-1000 µg/mL in Acetonitrile(±)(*)

Acetamiprid [100 µg/mL]	Aldicarb [1000 µg/mL]
Buprofezin [20 µg/mL]	Carbaryl [50 µg/mL]
Cyprodinil [250 µg/mL]	Dimethomorph [50 µg/mL]
Etiozazole [20 µg/mL]	Flonicamid [50 µg/mL]
Imidacloprid [20 µg/mL]	Iprodione [1000 µg/mL]
Mevinphos [50 µg/mL]	Novaluron [50 µg/mL]
Azoxystrobin [20 µg/mL]	Boscalid [20 µg/mL]
Carbofuran [20 µg/mL]	Chlorantraniliprole [20 µg/mL]
Dinotefuran [100 µg/mL]	Etofenprox [50 µg/mL]
Hexythiazox [10 µg/mL]	Imazalil [50 µg/mL]
Malathion [20 µg/mL]	Methiocarb [20 µg/mL]
Phosmet [20 µg/mL]	Piperonyl butoxide [200 µg/mL]

Canada Pesticide Mixture 2A

DRE-GA09001037AL 1ml
Canada Pesticide Mixture 2A 100 µg/mL in Acetonitrile(±)
DRE-GS09001038AL 5x1ml
Canada Pesticide Mixture 2A 100 µg/mL in Acetonitrile(±)

acephate	allethrin
buprofezin	chlorfenapyr
cypermethrin (mix of isomers)	cyprodinil
endosulfan II	endosulfan sulfate
etridiazole	fenoxycarb
fipronil	fludioxonil
kresoxim methyl	methoprene (mix of isomers)
phenothrin	phosdrinTM (mevinphos)
Propiconazol (mix of isomers)	pyraclostrobin
tetramethrin	trifloxystrobin
baythroid (mix of 4 isomers)	bifenthrin
chlorpyrifos	coumaphos
deltamethrin	endosulfan I
ethofenprox	etiozazole
fensulfothion	fenvalerate (mix of diastereoisomers)
iprodione	Kinoprene
MGK-264 - isomer a	permethrin (mix of isomers)
Pirimicarb	prallethrin
pyridaben	resmethrin

Canada Pesticide Mixture 3

DRE-GA09001042TO 1ml
Canada Pesticide Mixture 3 100 µg/mL in Toluene(±)(*)

acequinocyl	azadirachtin (Technical)
azoxystrobin	carbaryl
clofentezine	clothianidin
cyantraniliprole	dodemorph
fluopyram	hexythiazox
pentachloronitrobenzene	pyrethrin (mix of isomers)
tebufenozide	teflubenzuron

Canada Pesticide Mixture 4 ver. 2

DRE-A50000074AL 1ml
Canada Pesticide Mixture 4 ver. 2 20-500 µg/mL in Acetonitrile(±)
DRE-S50000074AL 5x1ml
Canada Pesticide Mixture 4 ver. 2 20-500 µg/mL in Acetonitrile(±)

Acequinocyl [30 µg/mL]	alpha-Endosulfan [200 µg/mL]
Benzovindiflupyr [20 µg/mL]	beta-Endosulfan [50 µg/mL]
Bifenazate [20 µg/mL]	Cyfluthrin [200 µg/mL]
Deltamethrin [500 µg/mL]	Endosulfan-sulfate [50 µg/mL]
Fenoxycarb [20 µg/mL]	Fenpyroximate (E/Z) [20 µg/mL]
Fenvalerate [100 µg/mL]	Permethrin [500 µg/mL]
Quintozene [20 µg/mL]	Thiophanate-methyl [50 µg/mL]

Canada Pesticide Mixture 5 ver. 2

DRE-A5000075EA 1ml
 Canada Pesticide Mixture 5 ver. 2 3000 µg/mL in Ethyl Acetate(±)
 DRE-S5000075EA 5x1ml
 Canada Pesticide Mixture 5 ver. 2 3000 µg/mL in Ethyl Acetate(±)

Oxamyl Spiromesifen

Carbamate Pesticide Mixture 673

DRE-A50000673ME 1ml
 Carbamate Pesticide Mixture 673 79 µg/mL (100 mg/kg) in Methanol(±)

aminocarb	methomyl
carbaryl	metolcarb
Pirimicarb	carbofuran
propoxur	fenobucarb
methiocarb	Indoxacarb (mixture of isomers)
aldicarb	aldicarb sulfone
aldicarb sulfoxide	isoprocarb
3-hydroxycarbofuran	

Carbamate Pesticides Mixture 638

DRE-A50000638ME 1ml
 Carbamate Pesticides Mixture 638 100 µg/mL in Methanol(±)

aldicarb	aldicarb sulfone
aldicarb sulfoxide	carbaryl
carbofuran	3-hydroxycarbofuran
methiocarb	methomyl
1-naphthol	oxamyl
propoxur	

Chinese Pharmacopoeia Pesticides Mixture 144

DRE-A50000144AL 1ml
 Chinese Pharmacopoeia Pesticides Mixture 144 100 µg/mL in Acetonitrile(±)

Aldrin	Chlordimeform free base
4,4'-DDE	2,4'-DDT
Dicofol	Dieldrin
Endosulfan-sulfate	Ethoprophos
Fipronil Sulfide	Fipronil Sulfone
beta-HCH	delta-HCH
Isofenphos-methyl	Monocrotophos
Parathion-methyl	Phorate
Terbufos	4,4'-DDD
Coumaphos	Demeton (O+S)
4,4'-DDT	beta-Endosulfan
alpha-Endosulfan	Fipronil
Fenamiphos	alpha-HCH
Fipronil-desulfinyl	Isocarbofos
gamma-HCH	Parathion-ethyl
Nitrofen	Sulfotep
Phospholan-methyl	

Chinese Pharmacopoeia Pesticides Mixture 146

DRE-A50000146AL 1ml
 Chinese Pharmacopoeia Pesticides Mixture 146 100 µg/mL in Acetonitrile(±)(*)

Aldicarb	Aldicarb-sulfone
Carbofuran	Carbofuran-3-hydroxy
Coumaphos	Demeton (O+S)
Fenamiphos	Fenamiphos-sulfone
Isazofos	Isocarbofos
Metsulfuron-methyl	Monocrotophos
Phorate-sulfoxide	Phosfolan
Terbufos-sulfone	Terbufos-sulfoxide
Aldicarb-sulfoxide	Cadusafos
Chlordimeform free base	Chlorsulfuron
Ethametsulfuron-methyl	Ethoprophos
Fenamiphos-sulfoxide	Fonofos
Isofenphos-methyl	Methamidophos
Phorate	Phorate-sulfone
Phosphamidon	Sulfotep

Canada Pesticide Mixture 6 ver. 2

DRE-A5000076IT 1ml
 Canada Pesticide Mixture 6 ver. 2 500-2000 µg/mL in Toluene:Isooctane(±)
 DRE-S5000076IT 5x1ml
 Canada Pesticide Mixture 6 ver. 2 500-2000 µg/mL in Toluene:Isooctane(±)

Kinoprene [500 µg/mL] Methoprene [2000 µg/mL]

Carbamate Pesticides Mixture 154 for HJ 827-2017

DRE-A50000154ME 1ml
 HJ 827-2017 Carbamate Pesticides Mixture 154 50-1000 µg/mL in Methanol(±)

Bendiocarb [200 µg/mL]	Carbofuran [50 µg/mL]
2,3,5-Trimethacarb [200 µg/mL]	Fenobucarb [100 µg/mL]
Propoxur [200 µg/mL]	Isoprocab [50 µg/mL]
Methiocarb [50 µg/mL]	Carbofuran-3-hydroxy [200 µg/mL]
Promecarb [100 µg/mL]	Metolcarb [100 µg/mL]
Pirimicarb [50 µg/mL]	Methomyl [200 µg/mL]
Carbaryl [200 µg/mL]	Methomyl-oxime [1000 µg/mL]

Chinese Pharmacopoeia Organochlorine Pesticides Mixture 114

DRE-A50000114IT 1ml
 Chinese Pharmacopoeia Organochlorine Pesticides Mixture 114 100-200 µg/mL in Isooctane:Toluene (±)

Aldrin [100 µg/mL]	beta-Endosulfan [100 µg/mL]
alpha-HCH [100 µg/mL]	beta-HCH [200 µg/mL]
Quintozene [100 µg/mL]	Heptachlor [100 µg/mL]
4,4'-DDE [100 µg/mL]	4,4'-DDD [200 µg/mL]
alpha-Endosulfan [100 µg/mL]	Endosulfan-sulfate [100 µg/mL]
delta-HCH [100 µg/mL]	Hexachlorobenzene [100 µg/mL]
2,4'-DDT [200 µg/mL]	gamma-HCH [100 µg/mL]
Endrin [200 µg/mL]	4,4'-DDT [100 µg/mL]
oxy-Chlordane [100 µg/mL]	Dieldrin [100 µg/mL]

Chinese Pharmacopoeia Pesticides Mixture 145

DRE-A50000145AL 1ml
 Chinese Pharmacopoeia Pesticides Mixture 145 40-100 µg/mL in Acetonitrile(±)(*)

Aldrin [100 µg/mL]	Chlordimeform free base [40 µg/mL]
4,4'-DDE [100 µg/mL]	2,4'-DDT [100 µg/mL]
Dicofol [100 µg/mL]	Dieldrin [100 µg/mL]
Endosulfan-sulfate [100 µg/mL]	Ethoprophos [40 µg/mL]
Fipronil Sulfide [40 µg/mL]	Fipronil Sulfone [40 µg/mL]
beta-HCH [100 µg/mL]	delta-HCH [100 µg/mL]
Isofenphos-methyl [40 µg/mL]	Monocrotophos [60 µg/mL]
Parathion-methyl [40 µg/mL]	Phorate [40 µg/mL]
Terbufos [40 µg/mL]	4,4'-DDD [100 µg/mL]
Coumaphos [100 µg/mL]	Demeton (O+S) [40 µg/mL]
4,4'-DDT [100 µg/mL]	beta-Endosulfan [100 µg/mL]
alpha-Endosulfan [100 µg/mL]	Fipronil [40 µg/mL]
Fenamiphos [40 µg/mL]	alpha-HCH [100 µg/mL]
Fipronil-desulfinyl [40 µg/mL]	Isocarbofos [100 µg/mL]
gamma-HCH [100 µg/mL]	Parathion-ethyl [40 µg/mL]
Nitrofen [100 µg/mL]	Sulfotep [40 µg/mL]
Phospholan-methyl [60 µg/mL]	

Chinese Pharmacopoeia Pesticides Mixture 147

DRE-A50000147AL 1ml
 Chinese Pharmacopoeia Pesticides Mixture 147 20-100 µg/mL in Acetonitrile(±)(*)

Aldicarb [100 µg/mL]	Aldicarb-sulfone [100 µg/mL]
Carbofuran [100 µg/mL]	Carbofuran-3-hydroxy [100 µg/mL]
Coumaphos [100 µg/mL]	Demeton (O+S) [40 µg/mL]
Fenamiphos [40 µg/mL]	Fenamiphos-sulfone [40 µg/mL]
Isazofos [20 µg/mL]	Isocarbofos [100 µg/mL]
Metsulfuron-methyl [100 µg/mL]	Monocrotophos [60 µg/mL]
Phorate-sulfoxide [40 µg/mL]	Phosfolan [60 µg/mL]
Terbufos-sulfone [40 µg/mL]	Terbufos-sulfoxide [40 µg/mL]
Aldicarb-sulfoxide [100 µg/mL]	Cadusafos [40 µg/mL]
Ethametsulfuron-methyl [100 µg/mL]	Ethoprophos [40 µg/mL]
Isofenphos-methyl [40 µg/mL]	Methamidophos [100 µg/mL]
Phorate [40 µg/mL]	Phorate-sulfone [40 µg/mL]
Phosphamidon [100 µg/mL]	Sulfotep [40 µg/mL]
Chlordimeform free base [40 µg/mL]	Chlorsulfuron [100 µg/mL]
Fenamiphos-sulfoxide [40 µg/mL]	Fonofos [40 µg/mL]

Chinese Pharmacopoeia Pesticides Mixture 352

DRE-A50000352AL 1ml
Chinese pharmacopoeia Pesticides Mixture 352 100 µg/mL in Acetonitrile(±)

2,4'-DDT	4,4'-DDD
Aldrin	alpha-Endosulfan
beta-HCH	Chlordimeform free base
Demeton (O+S)	Dicofol
Ethoprophos	Fenamiphos
Fipronil Sulfone	Fipronil-desulfinyl
Isofenphos-methyl	Monocrotophos
Parathion-ethyl	Parathion-methyl
Sulfotep	Terbufos
4,4'-DDE	4,4'-DDT
alpha-HCH	beta-Endosulfan
Coumaphos	delta-HCH
Dieldrin	Endosulfan-sulfate
Fipronil	Fipronil Sulfide
gamma-HCH	Isocarbofos
Nitrofen	o,p'-Dicofol
Phorate	Phospholan-methyl

Chlorinated Pesticides Mixture

DRE-A50000274TH 1ml
Chlorinated Pesticides Mixture 2000 µg/mL in Toluene/Hexane(±)

Aldrin	a-BHC
d-BHC	a-Chlordane
4,4'-DDE	4,4'-DDT
Endosulfan II	Endosulfan sulfate
Endrin ketone	Heptachlor
b-BHC	g-BHC
g-Chlordane	4,4'-DDD
Dieldrin	Endosulfan I
Endrin	Endrin aldehyde
Heptachlor epoxide (Isomer B)	Methoxychlor

Colorado Residual Pesticide Mixture

DRE-S50000081AC 5x1ml
Colorado Residual Pesticide Mixture 100 µg/mL in Acetone(±)(*)

Strobane	Aldrin
Phosphamidon	Methamidophos
gamma-HCH	HCH (BHC) (technical)
4,4'-DDT	4,4'-DDD
Fenoprop	Pentachlorophenol
Nitrofen	Dinoseb
2-Methyl-4,6-dinitrophenol	Daminozide
Parathion-ethyl	Parathion-methyl
Chlorobenzilate	Mevinphos
Endrin	Dieldrin
2,4-D-iso-octyl ester (technical)	Leptophos
Binapacryl	Hexachlorobenzene
Pyrinuron	Heptachlor
1,2-Dibromo-3-chloropropane	2,4,5-Trichlorophenoxyacetic acid
Captafol	2,4,5-Trichlorophenol
4-Chloranil	Fluoroacetamide
2-Ethyl-1,3-hexandiol	Safrole
MGK 11	EPN
Monocrotophos	Schradan
Chlordimeform free base	Mirex
Chlordecone	

EPA Method 505 Organohalide Pesticide Mixture 388

DRE-A50000388ME 1ml
EPA Method 505 Organohalide Pesticide Mixture 388 200 µg/mL in Methanol(±)

Alachlor	Atrazine
Aldrin	cis-Chlordane (alpha-Chlordane)
Endrin	gamma-HCH (Lindane)
Hexachlorobenzene	Hexachlorocyclopentadiene
Simazine	Methoxychlor (DMTD)
trans-Chlordane (gamma-Chlordane)	Dieldrin
Heptachlor	Heptachlor-exo-epoxide
cis-Nonachlor	trans-Nonachlor

Chlorinated Acids Mix 1

DRE-XA05550100AL 1ml
Chlorinated Acids Mix 1 100 µg/mL in Acetonitrile

2,4-D	Acifluorfen
Bentazone	Chloramben
Dicamba	Dichlorprop
Fenoprop	Picloram

Chlorinated Pesticide Mix 1

DRE-YA05080100MB 1ml
Chlorinated Pesticide Mix 1 1000 µg/mL in Methyl-tert-butyl ether(±)

4,4'-DDD	4,4'-DDE
alpha-Endosulfan	alpha-HCH
delta-HCH	Dieldrin
Endrin-aldehyde	gamma-HCH
Methoxychlor	Aldrin
4,4'-DDT	beta-HCH
beta-Endosulfan	Endrin
Endosulfan-sulfate	Heptachlor-exo-epoxide (cis-isomer B)
Heptachlor	

Colorado Pesticide Mixture 260

DRE-GA09000260AL 1ml
Colorado Pesticide Mixture 260 100 µg/mL in Acetonitrile(±)(*)

DRE-GS09000260AL 5x1ml
Colorado Pesticide Mixture 260 100 µg/mL in Acetonitrile(±)(*)

abamectin	azoxystrobin
bifenazate	etoxazole
imazalil	imidacloprid
malathion	permethrin (mixture of isomers)
spinosad (Spinosyn A & D)	spiromesifen
spirotetramat	Sythane TM
tebuconazole (Folicur)	

Diquat and Paraquat Mix 1

DRE-YA05490100WA 1ml
Diquat and Paraquat Mix 1 1000 µg/mL in Water

Diquat dibromide hydrate	Paraquat dichloride hydrate
--------------------------	-----------------------------

DZ/T 0064.72-2021 Organophosphorus Pesticide Mixture 696

DRE-A50000696AC 1ml
DZ/T 0064.72-2021 Organophosphorus Pesticide Mixture 696
1000 µg/mL in Acetone(±)

Dichlorvos	Phorate
Dimethoate	Parathion-methyl
Malathion	Chlorpyrifos
Parathion-ethyl	

EN 12918 Organophosphorus Mixture 439/440

DRE-B50000439AC 10ml
EN 12918 Organophosphorus Mixture 439 10 µg/mL in Acetone(±)

DRE-A50000440AC 1ml
EN 12918 Organophosphorus Mixture 440 100 µg/mL in Acetone(±)

Azinphos-ethyl	Azinphos-methyl
Chlorpyrifos (Chlorpyrifos-ethyl)	Chlorpyrifos methyl
Dimethoate	Fenitrothion
Mevinphos	Parathion (Parathion-ethyl)
Propetamphos	Triadimefon
Bromophos (Bromophos-methyl)	Chlorfenvinphos
Diazinon	Dichlorvos
Fenthion	Malathion
Parathion-methyl	Phosalone
Triazophos	

EPA Method 507 Pesticide Mixture 1

DRE-A50000461MB 1ml
 EPA Method 507 Pesticide Mixture 1 1000 µg/mL in Methyl tert Butyl Ether(‡)

Atrazine	Diphenamid
EPTC	Ethoprophos
Mevinphos	Prometryn
Propazine	Terbutryn
Triadimefon	

EPA Method 507 Pesticide Mixture 3

DRE-A50000463MB 1ml
 EPA Method 507 Pesticide Mixture 3 1000 µg/mL in Methyl tert Butyl Ether(‡)

Propachlor	Trifluralin
Benfluralin	Profluralin
Isopropalin	Pendimethalin
Oxadiazon	Oxyfluorfen

EPA Method 515 Herbicide Mixture

DRE-YS09000050AC 5x1ml
 EPA Method 515 Herbicide Mixture 100-1000 µg/mL in Acetone(‡)(*)

MCPP acid [10000 µg/mL]	MCPA acid [10000 µg/mL]
pentachlorophenol [100 µg/mL]	2,4,5-TP (Silvex) [100 µg/mL]
tetrachloroterephthalic acid [100 µg/mL]	4-nitrophenol [100 µg/mL]
2,4-D [100 µg/mL]	acifluorfen [100 µg/mL]
bentazon [100 µg/mL]	chloramben [100 µg/mL]
3,5-dichlorobenzoic acid [100 µg/mL]	Dicamba [100 µg/mL]
2,4,5-T [100 µg/mL]	picloram [100 µg/mL]
dichlorprop (2,4-DP) [100 µg/mL]	dinoseb [100 µg/mL]
dalapon [100 µg/mL]	2,4-DB [100 µg/mL]

EPA Method 515.2 Herbicide Mixture 458

DRE-A50000458ME 1ml
 EPA Method 515.2 Herbicide Mixture 458 100-1000 µg/mL in Methanol(‡)

Acifluorfen methyl ester [200 µg/mL]	Bentazon methyl [1000 µg/mL]
2,4-D methyl ester [100 µg/mL]	2,4-DB methyl ester [1000 µg/mL]
Dicamba-methyl ester [300 µg/mL]	Picloram methyl ester [300 µg/mL]
Fenoprop-methyl ester [100 µg/mL]	

EPA Method 515.4 Herbicide Mixture 409

DRE-A50000409MB 1ml
 EPA Method 515.4 Herbicide Mixture 409
 10-100 µg/mL in Methyl tert Butyl Ether(‡)

Acifluorfen methyl ester [50 µg/mL]	Bentazon methyl [100 µg/mL]
Dalapon methyl ester [100 µg/mL]	2,4-DB methyl ester [100 µg/mL]
Methyl-3,5-dichlorobenzoate [50 µg/mL]	Dichlorprop methyl ester [100 µg/mL]
Picloram methyl ester [50 µg/mL]	2,4,5-T methyl ester [25 µg/mL]
Chloramben methyl ester [50 µg/mL]	2,4-D methyl ester [100 µg/mL]
Dacthal [100 µg/mL]	Dicamba-methyl ester [50 µg/mL]
Dinoseb methyl ether [100 µg/mL]	Pentachloroanisole [10 µg/mL]
Fenoprop-methyl ester [25 µg/mL]	Quinlorac methyl ester [50 µg/mL]

EPA Method 525.3 Organophosphate Pesticide Mixture

DRE-GS09000342AC 5x1ml
 EPA Method 525.3 Organophosphate Pesticide Mixture
 500 µg/mL in Acetone(‡)

chlorfenvinphos (E/Z mixture)	chlorpyrifos
dimethipin	disulfoton
parathion	methyl parathion
phosphamidon	profenofos
dichlorvos	diisopropyl methylphosphonate
ethion	ethoprophos (prophos)
phosdrinTM (mevinphos)	phorate
tetrachlorvinphos (ISO)	tribufos

EPA Method 507 Pesticide Mixture 2

DRE-A50000462MB 1ml
 EPA Method 507 Pesticide Mixture 2 1000 µg/mL in Methyl tert Butyl Ether(‡)

Alachlor	Atraton
Bromacil	Butylate
Chlorpropham	Hexazinone
Molinate	Propyzamide (Pronamide)
Tetrachlorvinphos	Tricyclazole

EPA Method 507 Pesticide Mixture Kit 465

DRE-K50000465MB 1ea
 EPA Method 507 Pesticide Mixture Kit 465 1000 µg/mL in Methyl tert Butyl Ether(‡)

DRE-A50000461MB	1x1ml
EPA Method 507 Pesticide Mixture 1	1000 µg/mL in MtBE
DRE-A50000462MB	1x1ml
EPA Method 507 Pesticide Mixture 2	1000 µg/mL in MtBE
DRE-A50000463MB	1x1ml
EPA Method 507 Pesticide Mixture 3	1000 µg/mL in MtBE

EPA Method 515.2 Herbicide Mixture 402

DRE-A50000402ME 1ml
 EPA Method 515.2 Herbicide Mixture 402 100-1000 µg/mL in Methanol(‡)(*)

Acifluorfen [200 µg/mL]	Bentazon [1000 µg/mL]
2,4-D [100 µg/mL]	2,4-DB [1000 µg/mL]
Dicamba [300 µg/mL]	Picloram [300 µg/mL]
Fenoprop (Silvex) [100 µg/mL]	

EPA Method 515.2 Underivatized Mixture 404

DRE-A50000404ME 1ml
 EPA Method 515.2 Underivatized Mixture 404 100-500 µg/mL in Methanol(‡)

DCPA Diacid [100 µg/mL]	3,5-Dichlorobenzoic acid [500 µg/mL]
Dichlorprop [100 µg/mL]	Dinoseb [200 µg/mL]
Pentachlorophenol [100 µg/mL]	2,4,5-T [100 µg/mL]

EPA Method 525.3 Organonitrogen Pesticide Mixture (A-M)

DRE-S50000480AC 5x1ml
 EPA Method 525.3 Organonitrogen Pesticide Mixture 1
 500 µg/mL in Acetone(‡)

2,4-Dinitrotoluene	2,6-Dinitrotoluene
Atraton	Atrazine
Butylate	Butylhydroxytoluene
Cycloate	Diethyltoluamide (DEET)
Etridiazole	Fenarimol
Metolachlor	MGK 264
Alachlor	Ametryn
Bromacil	Butachlor
Chlorpropham	Cyanazine
Diphenamid	EPTC
Fluridone	Hexazinone

EPA Method 528 Phenol Calibration Mixture 389

DRE-A50000389DI 1ml
 EPA Method 528 Phenol Calibration Mixture 389
 2000 µg/mL in Dichloromethane(‡)

4-Chloro-3-methylphenol	2-Chlorophenol
2,4-Dichlorophenol	2,4-Dimethylphenol
2-Methyl-4,6-dinitrophenol	2,4-Dinitrophenol
2-Methylphenol	2-Nitrophenol
4-Nitrophenol	Pentachlorophenol
Phenol	2,4,6-Trichlorophenol

EPA Method 531.1 Carbamate Pesticide Mixture

DRE-GA09000948ME 1ml
 EPA Method 531.1 Carbamate Pesticide Mixture 100 µg/mL in Methanol(±)(*)

aldicarb	aldicarb sulfone
aldicarb sulfoxide	carbaryl
carbofuran	methiocarb
methomyl	oxamyl
propoxur	3-hydroxycarbofuran

EPA Method 622.1 Pesticide Mixture 392

DRE-A50000392MB 1ml
 EPA Method 622.1 Pesticide Mixture 392
 1000 µg/mL in Methyl tert Butyl Ether(±)

Aspon	Dichlofenthion
Famphur	Fenitrothion
Fonofos	Phosmet
Thionazin	

EPA Method 515.4 Herbicide Mixture 408

DRE-A50000408AC 1ml
 EPA Method 515.4 Herbicide Mixture 408 10-100 µg/mL in Acetone(±)(*)

Acifluorfen [50 µg/mL]	Bentazon [100 µg/mL]
Dalapon [100 µg/mL]	2,4-DB [100 µg/mL]
Dicamba [50 µg/mL]	3,5-Dichlorobenzoic acid [50 µg/mL]
Pentachlorophenol [10 µg/mL]	Picloram [50 µg/mL]
Quinclorac [50 µg/mL]	2,4-D [100 µg/mL]
Chloramben [50 µg/mL]	DCPA monoacid [50 µg/mL]
DCPA Diacid [50 µg/mL]	Dinoseb [100 µg/mL]
Dichlorprop [100 µg/mL]	Fenoprop (Silvex) [25 µg/mL]
2,4,5-T [25 µg/mL]	

EPA Method 614 Organophosphorus Pesticides Mixture

DRE-A50000275AH 1ml
 EPA Method 614 Organophosphorus Pesticides Mixture
 1000 µg/mL in Acetone/Hexane(±)

Azinphos-methyl	Demeton (mixed isomers)
Diazinon	Disulfoton
Ethion	Malathion
Parathion-ethyl	Parathion-methyl

EPA Method 8080A Organochlorine Pesticide QC Standard Mixture 467

DRE-A50000467AC 1ml
 EPA Method 8080A Organochlorine Pesticide QC Standard Mixture 467
 20-100 µg/mL in Acetone(±)

Aldrin [20 µg/mL]	alpha-HCH [20 µg/mL]
gamma-HCH (Lindane) [20 µg/mL]	4,4'-DDD (TDE) [100 µg/mL]
Dieldrin [20 µg/mL]	Endosulfan-alpha [20 µg/mL]
Endrin [100 µg/mL]	Endrin aldehyde [20 µg/mL]
Methoxychlor (DMTD) [20 µg/mL]	delta-HCH [20 µg/mL]
beta-HCH [20 µg/mL]	4,4'-DDT [100 µg/mL]
4,4'-DDE [20 µg/mL]	Endosulfan-total sulfate [100 µg/mL]
Endosulfan-beta [100 µg/mL]	Heptachlor-exo-epoxide [20 µg/mL]
Heptachlor [20 µg/mL]	

EPA Method 8141 Organophosphorus Pesticide Mixture 427

DRE-A50000427AH 1ml
 EPA Method 8141 Organophosphorus Pesticide Mixture 427
 200 µg/mL in n-Hexane:Acetone(±)

Dimethoate	EPN
Malathion	Monocrotophos
O,O-TEPP	Parathion-ethyl
Sulfotep	

EPA Method 608 Organochlorine Pesticide Mixture 391

DRE-A50000391IO 1ml
 EPA Method 608 Organochlorine Pesticide Mixture 391
 20 µg/mL in Isooctane(±)

Aldrin	alpha-HCH
delta-HCH	4,4'-DDD (TDE)
Dieldrin	Endosulfan-alpha
Endrin	Endrin aldehyde
beta-HCH	gamma-HCH (Lindane)
4,4'-DDE	4,4'-DDT
Endosulfan-beta	Endosulfan-total sulfate
Heptachlor	Heptachlor-exo-epoxide

EPA Method 1311 TCLP Methylated Herbicide Spiking Mixture 400

DRE-A50000400ME 1ml
 EPA Method 1311 TCLP Methylated Herbicide Spiking Mixture 400
 2000 µg/mL in Methanol(±)

2,4-D methyl ester	Fenoprop-methyl ester
--------------------	-----------------------

EPA Method 525.2 Organochlorine Pesticides Mixture

DRE-A50000278AC 1ml
 EPA Method 525.2 Organochlorine Pesticides Mixture
 100 µg/mL in Acetone(±)

Alachlor	Aldrin
b-BHC	d-BHC
Chlorothalonil	Chloroneb
p,p'-DDE	p,p'-DDT
Endosulfan II	Endosulfan sulfate
Etridiazole	a-Chlordane
Heptachlor epoxide (Isomer B)	Methoxychlor
Simazine	trans-Nonachlor
Atrazine	a-BHC
g-BHC	Chlorobenzilate
Dacthal	p,p'-DDD
Dieldrin	Endosulfan I
Endrin	Endrin aldehyde
g-Chlordane	Heptachlor
cis-Permethrin	trans-Permethrin

EPA Method 8080A Organochlorine Pesticide Mixture 613

DRE-A50000613TH 1ml
 EPA Method 8080A Organochlorine Pesticide Mixture 613
 1000 µg/mL in Hexane:Toluene(±)

o,p'-DDD	o,p'-DDT
o,p'-DDE	

EPA Method 8080A/8081 Organochlorine Pesticide Mixture 616

DRE-A50000616TH 1ml
 EPA Method 8080A/8081 Organochlorine Pesticide Mixture 616
 200 µg/mL in Hexane:Toluene(±)

hexachlorobenzene	a-BHC
b-BHC	d-BHC
g-BHC	aldrin
heptachlor	heptachlor epoxide isomer B
p,p'-DDD	p,p'-DDE
p,p'-DDT	dieldrin
endrin	o,p'-DDT

EPA Method 8141 Organophosphorus Pesticide Mixture 428

DRE-A50000428HE 1ml
 EPA Method 8141 Organophosphorus Pesticide
 Mixture 428 200 µg/mL in n-Hexane(±)

Carbophenothion	Chlorfenvinphos
Dioxathion	Ethion
Famphur	Azinphos-ethyl
Leptophos	Phosmet
Phosphamidon	Terbufos

EPA Method 8141 Organophosphorus Pesticide Mixture 429

DRE-A50000429HE 1ml
EPA Method 8141 Organophosphorus Pesticide Mixture 429
200 µg/mL in n-Hexane(±)

Aspon	Chlorpyrifos methyl
Crotoxyphos	Dichlofenthion
Dicrotophos	Fenitrothion
Fonofos	Thionazin
Trichlorfon	

Fenthion D6 & Atrazine D5 Mixture 579

DRE-A50000579AC 1ml
Fenthion D6 & Atrazine D5 Mixture 579 50 µg/mL in Acetone(±)

fenthion-d6	atrazine-d5
-------------	-------------

GB 23200.113-2018 Group B 105 Pesticides

DRE-A50000093EA 1.5ml
GB 23200.113-2018 Group B 105 Pesticides 10 µg/mL in Ethyl acetate(±)

aldrin as chlorine	acrinathrin [ISO]
atrazine	baythroid (mixture four of isomers)
Benfluralin (Benefin)	bifenox
butachlor	butamifos
chlorfenvinphos (E/Z-mixture)	chloroneb
chlorpropham	chlorpyrifos
cyproconazole (diastereomers)	cyprodinil
diazinon	dibrom
dieldrin	Difenoconazole (isomeric mixt.)
dipropetryn	ethiolat
etoxazole	etridiazole
fenbuconazole	fenchlorphos
fipronil	fluazifop-butyl
Fluorodifen	fluquinconazole
iprodione	isazophos
isoprothiolane	lambda cyhalothrin
malathion	mefenacet
methoxychlor	methyl parathion
ametryne	atraton
beflubutamid	benalaxyl
biphenyl	Bromophos ethyl
carbofuran	chlorfenson
chlorobenzilate	chlorpyrifos-methyl
Command (clomazone)	coumaphos
danitol	desmetryn
diclofop-methyl	dicrotophos
diniconazole (E isomer)	diphenylamine
ethion	ethofumesate
etrimfos	famphur
fenitrothion	fenobucarb
flucythrinate	fludioxonil
Guthion Ethyl	Hexaconazole
isocarbophos	isofenphos-oxon
leptophos	malaoxon
methidathion	methoprene (mixture of isomers)
monolinuron	napropamide
Nitrofen	omethoate
pendimethalin	pentachloroaniline
phosfolan	phosmet
procymidone	profenofos
Propiconazol (mixture of isomers)	propyzamide (pronamide)
simazine	Systhane TM
terbutylazine	terbutryne
thionazine (zinphos)	Tokuthion®
trichloronate	paclobutrazol (isomeric mixture)
oxadixyl	phosalone
pentachloronitrobenzene	Pirimiphos-ethyl
phosphamidon	Propanil
prometryn	pyrimethanil
pyridaphenthion	tecnazene
tau-fluvalinate	tetraconazole
tetrachlorvinphos (Rabon)	trans-chlordane
tolclofos-methyl	

EPA Method 8270 BN Mixture 207

DRE-GS09000207DI 5x1ml
EPA Method 8270 BN Mixture 207 2000 µg/mL in Dichloromethane(±)

2-chloronaphthalene	1,2-dichlorobenzene
1,3-dichlorobenzene	1,4-dichlorobenzene
hexachlorobenzene	hexachlorobutadiene
hexachlorocyclopentadiene	hexachloroethane
1,2,4-trichlorobenzene	2,4-dinitrotoluene
2,6-dinitrotoluene	isophorone
nitrobenzene	azobenzene

GB 23200.100-2016 Pyrethroide Pesticide Mixture 677

DRE-A50000677TH 1ml
GB 23200.100-2016 Pyrethroide Pesticide Mixture 677
100 µg/mL in Toluene:Hexane(±)

bifenthrin	danitol
lambda cyhalothrin	permethrin (mixture of isomers)
baythroid (mixture four of isomers)	cypermethrin (mix of isomers)
tau-fluvalinate	fenvalerate (mixture of diastereoisomers)
deltamethrin	

GB 23200.121-2021 Pesticide Mixture 1

DRE-A50000721AL 1ml
GB 23200.121-2021 Pesticide Mixture 1 50 µg/mL in Acetonitrile(±)(*)

Aldicarb	Aldicarb-sulfone
Bendiocarb	Carbaryl
Chlorantranilprole	Chlorpropham
Fenobucarb	Fenothiocarb
Isoprocarb	Methiocarb
Methomyl	Metolcarb
Oxamyl	Phenmedipham
Pirimicarb-desmethyl	Pirimicarb-desmethyl-formamido
Propoxur	Prosulfocarb
Triallate	Benazolin-ethyl ester
Aldicarb-sulfoxide	Carbofuran-3-hydroxy
Carbofuran	Famoxadone
Diethofencarb	Iprovalicarb
Fenoxycarb	Methiocarb sulfoxide
Methiocarb sulfone	Oxadixyl
Oxadiargyl	Pirimicarb
Phosalone	Propamocarb free base
Promecarb	Pyridaphenthion
Pyraclostrobin	

GB 23200.121-2021 Pesticide Mixture 2

DRE-A50000722AL 1ml
GB 23200.121-2021 Pesticide Mixture 2 50 µg/mL in Acetonitrile(±)(*)

Amidosulfuron	Bensulfuron-methyl
Ethoxysulfuron	Etrifos
Halosulfuron-methyl	Iodosulfuron-methyl sodium
Metsulfuron-methyl	Propyrisulfuron
Triasulfuron	Tribenuron-methyl (technical)
Chlorsulfuron	Cinosulfuron
Fenpyrazamine	Flucetosulfuron
Mesosulfuron-methyl	Metazosulfuron
Pyrazosulfuron-ethyl	Thifensulfuron-methyl
Triflusulfuron-methyl	Tritosulfuron

GB 23200.121-2021 Pesticide Mixture 3

DRE-A50000723AL 1ml
GB 23200.121-2021 Pesticide Mixture 3 50 µg/mL in Acetonitrile(±)

Acephate	Anilofos
Chlorpyrifos-methyl	Coumaphos
Demeton-S-methyl	Demeton-S-methyl sulfone
Dimethoate	Disulfoton
Disulfoton-sulfoxide	Edifenphos
Ethoprophos	Fensulfothion-sulfone
Fenthion-sulfoxide	Formothion
Isocarbafos	Isofenphos-methyl
Methacrifos	Methidathion
Phenthoate	Phorate
Phosmet	Phosmetoxon
Profenofos	Quinalphos
Terbufos-sulfone	Terbufos-sulfoxide
Triazophos	Vamidathion

GB 23200.121-2021 Pesticide Mixture 4

DRE-A50000724AL 1ml
 GB 23200.121-2021 Pesticide Mixture 4 50 µg/mL in Acetonitrile(±)

Chlorfluazuron	Chlorobenzuron
Fenamiphos	Fenamiphos-sulfone
Fenamiphos-sulfoxide	Hexaflumuron
Imazalil	Novaluron
Phosfolan	Phospholan-methyl
Procymidone	Pyribenzoxim
Teflubenzuron	Triflumuron

Cadusafos	Chlorpyrifos
Demeton (O+S)	Demeton-S Sulfoxide
Demeton-S-methyl sulfoxide	Diazinon
Disulfoton-oxon-sulfon	Disulfoton-sulfone
EPN	Ethion
Fenthion	Fenthion-sulfone
Iprobenfos	Isazofos
Malaoxon	Malathion
Omethoate	Parathion-ethyl
Phorate-sulfone	Phorate-sulfoxide
Phoxim	Pirimiphos-methyl
Sulfotep	Terbufos
Thidiazuron	Tolclofos-methyl

GB 23200.121-2021 Pesticide Mixture 5

DRE-A50000725AL 1ml
 GB 23200.121-2021 Pesticide Mixture 5 50 µg/mL in Acetonitrile(±)(*)

Acetochlor	Alachlor
Benalaxyl	Benzovindiflupyr
Butachlor	Carboxin
Clethodim Sulfone	Clomazone
Cyazofamid-dessulfonamide	Cyflufenamid
Diffufenican	Dimethenamid
Ethiprole	Ethirimol
Fenhexamid	Fenoxanil
Fensulfothion-oxon-sulfone	Fipronil
Florasulam	Flubendiamide
Flumorph	Fluopicolide
Amisulbrom	Azinphos-methyl
Benzoximate	Boscalid
Chloridazon	Chromafenozide
Cyantraniliprole	Cyazofamid
Dichlorvos	Diclotophos
Dimethomorph	Dimoxystrobin
Fenamidone	Fenaminostrobin
Fensulfothion	Fensulfothion-oxon
Fipronil Sulfone	Flonicamid
Flufenacet	Flumetsulam
Fluopyram	

GB 23200.121-2021 Pesticide Mixture 6

DRE-A50000726AL 1ml
 GB 23200.121-2021 Pesticide Mixture 6 50 µg/mL in Acetonitrile(±)

(±)-Metamifop	Fluthiacet-methyl
Fonofos	Heptenophos
Isopyrazam	Isoxaflutole
Mepronil	Metalaxyl
Methoxyfenozide	Mevinphos
Oxamyl-oxime	Oxaziclomefone
Penthiopyrad	Phosphamidon
Probenazole	Propachlor
Propyzamide	Proquinazid
Saflufenacil	Sedaxane
Spirotetramat	Spirotetramat-enol-glucoside
Sulfoxaflor	Tebufenozide
Flutolanil	Fluxapyroxad
Hexazinone	Iprodione
Mandipropamid	Mefenacet
Metamitron	Methamidophos
Monocrotophos	Napropamide
Penflufen	Penoxsulam
Picolinafen	Pretilachlor
Propanil	Propisochlor
Pyridaben	Pyrimorph
Silthiofarm	S-Metolachlor
Spirotetramat-keto-hydroxy	Spirotetramat-mono-hydroxy
Tolfenpyrad	Zoxamide

GB 23200.121-2021 Pesticide Mixture 7

DRE-A50000727AL 1ml
 GB 23200.121-2021 Pesticide Mixture 7 50 µg/mL in Acetonitrile(±)

(E)-Nitenpyram	2-Diethylaminoethyl Hexanoate
Ametryn	Atrazine
Clothianidin	Cyanazine
Dinotefuran	Fipronil Sulfide
Imidaclothiz	Metribuzin
Pyraoxystrobin	Pyrimethanil
Spinosyn A	Spinosyn D
Thiabendazole	Thifluzamide
Acetamidprid	Ametoctradin
Bupirimate	Butralin
Cyprodinil	Diclobutrazol
Flurtamone	Imidacloprid
Pendimethalin	Prometryn
Simazine	Simetryn
tau-Fluvalinate	Terbuthylazine
Triflumizole-amino	

GB 23200.121-2021 Pesticide Mixture 8

DRE-A50000728MC 1ml
 GB 23200.121-2021 Pesticide Mixture 8 50 µg/mL in Acetonitrile:Methanol(±)(*)

(±)-Fenpropathrin	(E)-Fenpyroximate
Bioresmethrin	Bitertanol
Cyproconazole	Deltamethrin
Diniconazole	Dinocap
Epoxiconazole	Etofenprox
Fluazifop-butyl	Fluoroglycofen-ethyl
Kresoxim-methyl	Lactofen
Oxadiazon	Permethrin
Propaquizafop	Propargite
Pyriproxyfen	Pyriproxyfen
Spinetoram	Spirodiclofen
Thiacloprid	Thiamethoxam
Azoxystrobin	Bifenthrin
Coumoxystrobin	Cycloxydim
Diclofop methyl	Difenoconazole
Emamectin benzoate	Enoxastrobin
Etoazole	Fenvalerate
Isoprothiolane	Ivermectin
Methoprene	Myclobutanil
Picoxystrobin	Piperonyl butoxide
Pyraflufen-ethyl	Pyrethrin 1
Quinalofop-ethyl	Rotenone
Spiromesifen	Tebuconazole
Triadimenol	

GB 23200.121-2021 Pesticide Mixture 9

DRE-A50000729AL	1ml	
GB 23200.121-2021 Pesticide Mixture 9 50 µg/mL in Acetonitrile(‡)		
Avermectin B1a	BifenoX	
Chlorfenvinphos	Clethodim Sulfoxide	
Ethofumesate	Fenarimol	
Fenoxaprop-ethyl	Fenpropidin	
Fluazinam	Flucythrinate	
Flusilazole	Flutriafol	
Ipcnazole	Isoxaflutole-diketonitrile	
Metrafenone	Oxyfluorfen	
Phenamacril	Propiconazole	
Pyrisoxazole	Tetraconazole	
Tricyclazole	Trifloxystrobin	
Uniconazole	Chlordimeform free base	
Bromuconazole	Cyflumetofen	
Clofentezine	Fenbuconazole	
Fenazaquin	Fipronil-desulfinyl	
Fenpropimorph	Flumetralin	
Fludioxonil	Imibenconazole	
Hexaconazole	Metconazole	
Metazachlor	Penconazole	
Paclobutrazol	Pyridalyl	
Pyrametostrobin	Trichlorphon	
Triadimefon	Triticonazole	
Triflumizole		

GB 23200.121-2021 Pesticide Mixture Kit 1

DRE-K50000720	1ea	
GB 23200.121-2021 Pesticide Mixture Kit 1(‡)(*)		
DRE-A10065000AL-50	Albendazole 50 µg/mL in Acetonitrile	1x1ml
DRE-A10990000ME-50	Carbendazim 50 µg/mL in Methanol	1x1ml
DRE-A13497500AC-50	Fenothiocarb 50 µg/mL in Acetone	1x1ml
DRE-A15920000AL-50	Pencycuron 50 µg/mL in Acetonitrile	1x1ml
DRE-A16940000AL-50	Sethoxydim 50 µg/mL in Acetonitrile	1x1ml
DRE-A17605000AL-50	Tralkoxydim 50 µg/mL in Acetonitrile	1x1ml
DRE-A50000721AL	GB 23200.121-2021 Pesticide Mixture 1 50 µg/mL in Acetonitrile	1x1ml
DRE-A50000722AL	GB 23200.121-2021 Pesticide Mixture 2 50 µg/mL in Acetonitrile	1x1ml
DRE-A50000723AL	GB 23200.121-2021 Pesticide Mixture 3 50 µg/mL in Acetonitrile	1x1ml
DRE-A50000724AL	GB 23200.121-2021 Pesticide Mixture 4 50 µg/mL in Acetonitrile	1x1ml
DRE-A50000725AL	GB 23200.121-2021 Pesticide Mixture 5 50 µg/mL in Acetonitrile	1x1ml
DRE-A50000726AL	GB 23200.121-2021 Pesticide Mixture 6 50 µg/mL in Acetonitrile	1x1ml
DRE-A50000727AL	GB 23200.121-2021 Pesticide Mixture 7 50 µg/mL in Acetonitrile	1x1ml
DRE-A50000728MC	GB 23200.121-2021 Pesticide Mixture 8 50 µg/mL in Acetonitrile:Methanol	1x1ml
DRE-A50000729AL	GB 23200.121-2021 Pesticide Mixture 9 50 µg/mL in Acetonitrile	1x1ml
DRE-A50000730AL	GB 23200.121-2021 Pesticide Mixture 10 50 µg/mL in Acetonitrile	1x1ml

GB/T 14848-2017 Pesticides Mixture 522

DRE-A50000522TO	1ml	
GB/T 14848-2017 Pesticides Mixture 522 100 µg/mL in Toluene(‡)(*)		
Heptachlor	2,4-D ((2,4-Dichlorophenoxy)acetic Acid)	
Carbofuran	Aldicarb	
Dichlorvos	Parathion-methyl	
Malathion	Dimethoate	
Chlorpyrifos	Chlorothalonil	
Atrazine		

GB 23200.121-2021 Pesticide Mixture 10

DRE-A50000730AL	1ml	
GB 23200.121-2021 Pesticide Mixture 10 50 µg/mL in Acetonitrile(‡)(*)		
Buprofezin	Carfentrazone-ethyl	
Cyclosulfamuron	Cymoxanil	
Diuron	Flufenoxuron	
Hexythiazox	Indoxacarb	
Lufenuron	Metaflumizone	
Prochloraz desimidazole-amino	Prochloraz desimidazole-formylamino	
Thiophanate-methyl	Chlorotoluron	
Chlorimuron-ethyl	Dimepiperate	
Diflubenzuron	Fosthiazate	
Forchlorfenuron	Linuron	
Isoproturon	Prochloraz	
Molinate	Tebuthiuron	
Suffentrazone		

GB 31656.8-20 Organophosphorus Pesticide Mixture 694

DRE-A50000694AL	1ml	
GB 31656.8-20 Organophosphorus Pesticide Mixture 694 100-200 µg/mL in Acetonitrile(‡)		
Phoxim [200 µg/mL]	Propetamphos [100 µg/mL]	
Fenthion [200 µg/mL]	Malathion [100 µg/mL]	
Diazinon [100 µg/mL]	Trichlorfon [100 µg/mL]	
Dichlorvos [100 µg/mL]	Azamethiphos [100 µg/mL]	
Coumaphos [200 µg/mL]		

GB 31658.8-2021 Pyrethroid Mixture 695

DRE-A50000695AC	1ml	
GB 31658.8-2021 Pyrethroid Mixture 695 100 µg/mL in Acetone(‡)		
Deltamethrin	Bifenthrin	
Flucythrinate	tau-fluvalerate	
Tefluthrin	Fenvalerate	

GB 5009.168-2016 Fatty acid methyl esters

DRE-A50000712HP	5x1ml	
GB 5009.168-2016 Fatty acid methyl esters 200-400 µg/mL in n-Heptane(‡)(*)		
Arachidic acid methyl ester [200 µg/mL]	Methyl arachidonate [300 µg/mL]	
Capric acid methyl ester [200 µg/mL]	Caproic acid methyl ester [200 µg/mL]	
Me cis-10-pentadecenoate [200 µg/mL]	Me cis-11,14,17-eicosatrienoate [200]	
Me cis-8,11,14-eicosatrienoate [200]	Methyl dodecanoate [200 µg/mL]	
Y-Linolenic acid methyl ester [200 µg/mL]	Gondolic acid methyl ester [200 µg/mL]	
Linoelaidic methyl ester [200 µg/mL]	Linoleic acid methyl ester [400 µg/mL]	
Methyl docosahexaenoate [300 µg/mL]	Methyl myristate [200 µg/mL]	
Methyl stearate [400 µg/mL]	Myristoleic acid methyl ester [200 µg/mL]	
Methyl pentadecanoate [200 µg/mL]	Methyl tetracosanoate [200 µg/mL]	
Methyl undecanoate [200 µg/mL]	Butyric acid methyl ester [200 µg/mL]	
Behenic acid methyl ester [200 µg/mL]	Me cis-10-heptadecenoate [200 µg/mL]	
Caprylic acid methyl ester [200 µg/mL]	Me cis-13,16-docosadienoat [200µg/mL]	
Me cis-11,14-eicosadienoat [200µg/mL]	Erucic acid methyl ester [200 µg/mL]	
Methyl eicosapentaenoate [300 µg/mL]	Methyl heptadecanoate [200 µg/mL]	
Methyl heneicosanoate [200 µg/mL]	Methyl 9-octadecenoate [200 µg/mL]	
Linolenic acid methyl ester [300 µg/mL]	Methyl palmitoleate [200 µg/mL]	
Methyl palmitate [400 µg/mL]	Oleic acid methyl ester [400 µg/mL]	
Nervonic acid methyl ester [200 µg/mL]	Methyl tridecanoate [200 µg/mL]	
Methyl tricosanoate [200 µg/mL]		

GB/T 39665-2020 Pesticide Mixture 606

DRE-A50000606ME	1ml
GB/T 39665-2020 Pesticide Mixture 606 100 µg/mL in Methanol(±)(*)	
Carbaryl	Carbendazim
Diuron	Epoxiconazol
Flumioxazin	loxynil
Kresoxim-methyl	Linuron
Propargite	Propazine
Diallate	Dinoterb
Fenpropimorph	Fluazifop-P-butyl
Iprodione	Isoxaflutole
Molinate	Monocrotophos
Simazine	Tridemorph

GB/T 39665-2020 Pesticide Mixture 607

DRE-A50000607AC	1ml
GB/T 39665-2020 Pesticide Mixture 607 500-1000 µg/mL in Acetone(±)	
Alachlor [1000 µg/mL]	Aldrin [1000 µg/mL]
gamma-HCH [1000 µg/mL]	delta-HCH [1000 µg/mL]
trans-Chlordane [500 µg/mL]	Chlorothalonil [1000 µg/mL]
4,4'-DDD [1000 µg/mL]	2,4'-DDE [1000 µg/mL]
4,4'-DDT [1000 µg/mL]	1,4-Dichlorobenzene [1000 µg/mL]
Diphenylamine [1000 µg/mL]	α-Endosulfan [1000 µg/mL]
Endrin [1000 µg/mL]	Fenarimol [1000 µg/mL]
Heptachlor [1000 µg/mL]	Heptachlor-exo-epoxide [1000 µg/mL]
Pronamide [1000 µg/mL]	Vinclozolin [1000 µg/mL]
alpha-HCH [1000 µg/mL]	beta-HCH [1000 µg/mL]
Captafol [1000 µg/mL]	cis-Chlordane [500 µg/mL]
Chlozolinate [1000 µg/mL]	2,4'-DDD [1000 µg/mL]
4,4'-DDE [1000 µg/mL]	2,4'-DDT [1000 µg/mL]
Dicofol [1000 µg/mL]	Dieldrin [1000 µg/mL]
β-Endosulfan [1000 µg/mL]	Endosulfan-sulfate [1000 µg/mL]
Fenthion [1000 µg/mL]	Flusilazole [1000 µg/mL]
Hexachlorobenzene [1000 µg/mL]	Myclobutanil [1000 µg/mL]
Phosphamidon [1000 µg/mL]	

Group A 117 Pesticides for GB 23200.113-2018

DRE-A50000092EA 1.5ml
GB 23200.113-2018 Group A 117 Pesticides 10 µg/mL in Ethyl acetate(±)(*)

Tetramethrin	Dichlofenthion	Fenarimol	Bifenthrin	Dicofol
Allethrin	Propoxur	Isoprocarb	Tebupirifos	Pyriproxyfen
Dioxathion	(+)-trans-Permethrin	beta-Endosulfan	alpha-Endosulfan	Alachlor
Bromophos-methyl	Carbophenothion	Fenamidon	Mephosfolan	Ditalimfos
Bromfenvinphos	Deltamethrin	Pirimicarb	Bupirimate	Pyridaben
Methamidophos	Fenvalerate	Cypermethrin (technical)	Edifenphos	Metribuzin
Triadimenol	Triadimefon	Tebuconazole	Hexachlorobenzene	Piperonyl butoxide
alpha-HCH	beta-HCH	delta-HCH	gamma-HCH	Atrazine-desethyl
Tetradifon	Epoxiconazole	Imazalil	Penconazole	Phorate
Ethoprophos	Tribufos	2,4'-DDT	2,4'-DDE	Isoxathion
2,4'-DDD	4,4'-DDT	4,4'-DDE	4,4'-DDD	Fenthion
Dichlorvos	Dicloran	Dichlobenil	Monocrotophos	Iprobenfos
Oxyfluorfen	Aclonifen	Pretilachlor	Kresoxim-methyl	Phorate-sulfoxide
Metolachlor	Acetochlor	Boscalid	Mevinphos	Sulfotep
Formothion	Dimethoate	Piperophos	Tebufenpyrad	Fenthion-sulfoxide
Vinclozolin	Fosthiazate	Hexazinone	Propetamphos	EPN
Pirimiphos-methyl	Mepanipyrim	Quinoxifen	Endrin	Trifloxystrobin
Bromacil	Oxadiazon	Propazine	Molinate	Flutolanil
Fenthion-sulfone	Fensulfotiazon	Parathion-ethyl	Chlorthiophos	Cyflufenamid
Fonofos	Phorate-sulfone	Triazophos	Isofenphos-methyl	Bromopropylate
Methacrifos	Quinalphos	Paraoxon-ethyl	Thiobencarb	Triallate
Anilofos	Metalaxyl	Paraoxon-methyl	Terbufos-sulfone	Cycloate
Acephate	Profluralin	Pyrazophos	Fenothiocarb	Terbufos
Isofenphos	Ethalfuralin			

PESTICIDE MIXTURES

GB/T 5750.9-2006 Pyrethroids Pesticide Mixture 557

DRE-A50000557AL 1ml
GB/T 5750.9-2006 Pyrethroids Pesticide Mixture 557 50 µg/mL in Acetonitrile(±)

fenpropathrin	lambda cyhalothrin
permethrin (mixture of isomers)	cypermethrin (mix of isomers)
deltamethrin	fenvalerate (mixture of diastereoisomers)

Haloacetic Acid Mixture 540

DRE-A50000540MB 1ml
Haloacetic Acid Mixture 540 1000 µg/mL in Methyl tert Butyl Ether(±)

chloroacetic acid	dichloroacetic acid
trichloroacetic acid	

Halogenated Pesticide Mixture 562

DRE-A50000562PM 1ml
Halogenated Pesticide Mixture 562 100 µg/mL in Petroleum Ether(±)

o,p'-DDT	a-BHC
b-BHC	d-BHC
g-BHC	pentachloronitrobenzene
p,p'-DDT	p,p'-DDD
p,p'-DDE	

HJ 753-2015 Pyrethroid Pesticides Mixtures

DRE-A50000153AC 1ml
HJ 753-2015 Pyrethroid Pesticides Mixture 153 100 µg/mL in Acetone(±)
DRE-A50000611AC 1ml
HJ 753-2015 Pyrethroid Pesticides Mixture 611 1000 µg/mL in Acetone(±)

deltamethrin	fenvalerate (mixt. of diastereoisomers)
cypermethrin (mix of isomers)	lambda cyhalothrin
bifenthrin	tetramethrin
danitol	allethrin

HJ 1022-2019 Phenoxy-carboxylic Pesticides Mixture

DRE-A30000010AC 1ml
HJ 1022-2019 Phenoxy-carboxylic Pesticides Mixture 10
100 µg/mL in Acetone(±)

DRE-A30000014AC 1ml
HJ 1022-2019 Pesticides Mixture 14 500 µg/mL in Acetone(±)

2,4-DB	2,4-D
Dicamba	Dichlorprop
Fenoprop	MCPA
2,4,5-T	

HJ 1189-2021 Organophosphorus Pesticide Mixture 693

DRE-A50000693AC 1ml
HJ 1189-2021 Organophosphorus Pesticide Mixture 693
2000 µg/mL in Acetone(±)

Dichlorvos	Mevinphos
Sulfotep	Phorate
Fonofos	Probenfos
Chlorpyrifos-methyl	Phosphamidon
Fenitrothion	Malathion
Isofenphos-methyl	Isocarbafos
Fenamiphos	Triazophos
Demeton (O+S)	Ethoprophos
Terbufos	Diazinon
Dimethoate	Isazofos
Parathion-methyl	Chlorpyrifos
Parathion-ethyl	Bromophos-methyl
Phenthoate	Profenofos
Coumaphos	

HJ 350-2007 Organochlorine Pesticides Mixture 600

DRE-A50000600TH 1ml
HJ 350-2007 Organochlorine Pesticides Mixture 600
1000 µg/mL in Toluene:Hexane(±)

endrin ketone	heptachlor
heptachlor epoxide isomer B	hexachlorobenzene
hexachlorocyclopentadiene	isodrin
methoxychlor	toxaphene
chlordane (Mix of Isomers)	

HJ 350-2007 Organochlorine Pesticides Mixture 668

DRE-A50000668TH 1ml
HJ 350-2007 Organochlorine Pesticides Mixture 668
1000 µg/mL in Toluene:Hexane(±)

aldrin	a-BHC
d-BHC	cis-chlordane
p,p'-DDD	p,p'-DDE
dieldrin	endosulfan I
endrin	endrin aldehyde
b-BHC	g-BHC
trans-chlordane	1,2-dibromo-3-chloropropane
p,p'-DDT	di-allate (mixture of isomers)
endosulfan II	endosulfan sulfate
chlorobenzilate	

HJ 699-2014 Organochlorine Pesticides & Chlorobenzenes Mixture 11

DRE-A30000011AC 1ml
HJ 699-2014 Organochlorine Pesticides & Chlorobenzenes Mixture 11 100 µg/mL in Acetone

Aldrin	Dieldrin
Endrin-ketone	alpha-Endosulfan
Heptachlor	(±)-cis-Heptachlorepoxyde
trans-Chlordane	2,4'-DDD
4,4'-DDE	2,4'-DDT
Methoxychlor	alpha-HCH
delta-HCH	1,2,3-Trichlorobenzene
1,2,3,4-Tetrachlorobenzene	1,2,3,5-Tetrachlorobenzene
Hexachlorobenzene	Quintozone
Endrin	Endrin-aldehyde
beta-Endosulfan	Endosulfan-sulfate
(±)-trans-Heptachlorepoxyde	cis-Chlordane
4,4'-DDD	2,4'-DDE
4,4'-DDT	Dicofol
beta-HCH	gamma-HCH
1,2,4-Trichlorobenzene	1,3,5-Trichlorobenzene
1,2,4,5-Tetrachlorobenzene	Pentachlorobenzene

HJ 699-2014 Organochlorine Pesticides Mixture 517

DRE-A50000517HE 1ml
HJ 699-2014 Organochlorine Pesticides Mixture 517 10 µg/mL in n-Hexane(±)

alpha-HCH	beta-HCH
Quintozone	Heptachlor
(±)-trans-Heptachlorepoxyde	cis-Chlordane
beta-Endosulfan	Dieldrin
2,4'-DDE	Mitotane
4,4'-DDT	Methoxychlor
Endosulfan-sulfate	Dicofol
Lindane	delta-HCH
Aldrin	(±)-cis-Heptachlorepoxyde
trans-Chlordane	alpha-Endosulfan
Endrin	4,4'-DDE
4,4'-DDD	2,4'-DDT
Endrin-aldehyde	Endrin-ketone

HJ 961-2018, HJ 1026-2019 Internal Standards Mixture 498

DRE-A50000498ME 1ml
HJ 961-2018, HJ 1026-2019 Internal Standards Mixture 498
100 µg/mL in Methanol(±)(*)

Carbaryl D7 (naphthyl D7)	Methomyl D3
---------------------------	-------------

IRSA Method 5090 Pesticide Mixture

DRE-GA09000944AC 1ml
IRSA Method 5090 Pesticide Mixture 100 µg/mL in Acetone(±)

a-BHC	b-BHC
aldrin	heptachlor
p,p'-DDD	p,p'-DDE
endrin	atrazine
o,p'-DDE	isodrin
endosulfan I	endosulfan II
alachlor	g-BHC
d-BHC	trans-chlordane
cis-chlordane	dieldrin
p,p'-DDT	o,p'-DDT
o,p'-DDD	heptachlor epoxide isomer B
hexachlorobenzene	methoxychlor
pentachlorobenzene	

ISO 6468 Chlorinated Pesticide Mixture 362

DRE-A50000362IO 1ml
ISO 6468 Chlorinated Pesticide Mixture 362 10 µg/mL in Isooctane(‡)

alpha-HCH	beta-HCH
epsilon-HCH	2,4'-DDE
4,4'-DDD (TDE)	2,4'-DDT
Aldrin	Dieldrin
Heptachlor-endo-epoxide	Heptachlor-exo-epoxide
1,2,3-Trichlorobenzene	1,2,4-Trichlorobenzene
1,2,3,5-Tetrachlorobenzene	1,2,4,5-Tetrachlorobenzene
Quintozone (Pentachloronitrobenzene)	PCB 28
PCB 138	PCB 153
gamma-HCH (Lindane)	delta-HCH
4,4'-DDE	2,4'-DDD (o,p'-TDE)
4,4'-DDT	Methoxychlor (DMTD)
Endrin	Heptachlor
Endosulfan-alpha	Endosulfan-beta
1,2,5-Trichlorobenzene	1,2,3,4-Tetrachlorobenzene
Pentachlorobenzene	Hexachlorobenzene
PCB 52	PCB 101
PCB 180	PCB 194

ISO 10695:2000 Standard Mixture 366

DRE-B50000366AC 10ml
ISO 10695:2000 Standard Mixture 366 10 µg/mL in Acetone(‡)

Atrazine	Cyanazine
Metazachlor	Parathion (Parathion-ethyl)
Parathion-methyl	Pendimethalin
Propazine	Sebuthylazine
Simazine	Terbuthylazine
Trifluralin	Vinclozolin

Maryland Pesticide Mixture 2

DRE-A50000141AL 1ml
Maryland Pesticide Mixture 2 1000 µg/mL in Acetonitrile(‡)
DRE-S50000141AL 5x1ml
Maryland Pesticide Mixture 2 1000 µg/mL in Acetonitrile(‡)
DRE-A50000209AL 1ml
Maryland Pesticide Mixture 2 1000 µg/mL in Acetonitrile Second Source(‡)
DRE-S50000209AL 5x1ml
Maryland Pesticide Mixture 2 1000 µg/mL in Acetonitrile Second Source(‡)

Bifenazate	Bifenthrin
Cyfluthrin	Diazinon
Hexythiazox	Metalaxyl
Phosmet	Piperonyl butoxide
Boscalid	Chlorpyrifos
Fipronil	Flurprimidol
Paclbutrazol	Permethrin
Pyrethrins	Trifloxystrobin

Method DM 471 Pesticide Mixture

DRE-GA09000941AC 1ml
Method DM 471 Pesticide Mixture 100 µg/mL in Acetone(‡)

alachlor	aldrin
g-BHC	cis-chlordane
p,p'-DDT	p,p'-DDD
o,p'-DDT	dieldrin
a-BHC	b-BHC
trans-chlordane	o,p'-DDE
p,p'-DDE	o,p'-DDD
endrin	atrazine

ISO 10382:2002

PCB and Organochlorine Pesticide Mixture 369

DRE-A50000369IO 1ml

ISO 10382:2002 PCB and Organochlorine Pesticide Mixture 369
100 µg/mL in Isooctane(‡)

PCB 28	PCB 52
PCB 138	PCB 153
alpha-HCH	beta-HCH
Dieldrin	Endrin
Heptachlor-exo-epoxide	Endosulfan-alpha
2,4'-DDT	4,4'-DDD (TDE)
PCB 101	PCB 118
PCB 180	Hexachlorobenzene
gamma-HCH (Lindane)	Aldrin
Heptachlor	Heptachlor-endo-epoxide
4,4'-DDE	2,4'-DDD (o,p'-TDE)
2,4'-DDE	4,4'-DDT

Massachusetts Residual Pesticide Mixture

DRE-S50000048AL 5x1ml
Massachusetts Residual Pesticide Mixture
100 µg/mL in Acetonitrile(‡)

Imidacloprid	Imazalil
Sythane Tm	Bifenazate
Trifloxystrobin	Spiromesifen
Bifenthrin	Etoazole
Baythroid (mixture Four Of Isomers)	

Maryland Pesticide Mixture 1

DRE-A50000140AL 1ml
Maryland Pesticide Mixture 1 1000 µg/mL in Acetonitrile(‡)(*)
DRE-S50000140AL 5x1ml
Maryland Pesticide Mixture 1 1000 µg/mL in Acetonitrile(‡)(*)
DRE-S50000208AL 5x1ml
Maryland Pesticide Mixture 1 1000 µg/mL in Acetonitrile Second Source(‡)(*)

(E)-Fenpyroximate	Abamectin
Ancymidol	Azoxystrobin
Chlorantraniliprole	Dimethoate
Flonicamid	Fludioxonil
Myclobutanil	Propiconazole
Acetamiprid	Aldicarb
Carbaryl	Carbofuran
Ethephon	Etoazole
Imidacloprid	Methomyl
Thiacloprid	Thiamethoxam

Michigan Pesticide Mixture 2

DRE-A50000100AL 1ml
Michigan Pesticide Mixture 2 100 µg/mL in Acetonitrile(‡)(*)
DRE-S50000100AL 5x1ml
Michigan Pesticide Mixture 2 100 µg/mL in Acetonitrile(‡)(*)

Abamectin	Acetamiprid
Bifenthrin	Boscalid
Fenoxycarb	Fipronil
Imazalil	Imidacloprid
Permethrin	Prallethrin
Thiacloprid	Trifloxystrobin
Aldicarb	Azoxystrobin
Cyfluthrin	Cypermethrin (technical)
Flonicamid	Fludioxonil
Methiocarb	Myclobutanil
Pyrethrins	Spinosad

Michigan Residual Solvents Mixture 470

DRE-S50000470TN 5x1ml
Michigan Residual Solvents Mixture 470 100-1000 µg/mL in Triacetin(‡)

1,2-Dichloroethane [100 µg/mL]	2,2-Dimethylbutane [1000 µg/mL]
2-Methylpentane [1000 µg/mL]	3-Methylpentane [1000 µg/mL]
Benzene [100 µg/mL]	Chloroform [100 µg/mL]
Ethanol [1000 µg/mL]	Ethyl acetate [1000 µg/mL]
n-Heptane [1000 µg/mL]	n-Hexane [1000 µg/mL]
Trichloroethene [100 µg/mL]	Xylene (all isomers) [1000 µg/mL]
2,3-Dimethylbutane [1000 µg/mL]	2-Methylbutane [1000 µg/mL]
Acetone [1000 µg/mL]	Acetonitrile [1000 µg/mL]
Dichloromethane [1000 µg/mL]	Diethylether [1000 µg/mL]
Isopropyl alcohol [1000 µg/mL]	Methanol [1000 µg/mL]
n-Pentane [1000 µg/mL]	Toluene [1000 µg/mL]

Nevada Pesticide Mixture 62

DRE-GA09000062AL 1ml
Nevada Pesticide Mixture 62 100 µg/mL in Acetonitrile(‡)(*)
DRE-GS09000062AL 5x1ml
Nevada Pesticide Mixture 62 100 µg/mL in Acetonitrile(‡)(*)

acequinocyl	bifenazate
baythroid (mixture of isomers)	dimethomorph
pentachloronitrobenzene	spinosad (Spinosyn A & D)
cypermethrin (mix of isomers)	piperonyl butoxide
fenhexamid	flonicamid
fludioxonil	pyrethrin (mixture of isomers)
bifenthrin	captan
etoxazole	Systhane TM
thiamethoxam	trifloxystrobin
imidacloprid	abamectin
spinetoram (mixture of isomers)	spirotetramat

Nitrogen/Phosphorus Pesticide Mixture 925

DRE-GA090000925AC 1ml
Nitrogen/Phosphorus Pesticide Mixture 925 100 µg/mL in Acetone(‡)

alachlor	atrazine
propyzamide (pronamide)	propachlor
ametryne	atraton
prometon	tebuthiuron
butylate	chlorpropham
pebulate	terbacil
MGK-264 (mixture of isomers)	phosdrinTM (mevinphos)
dichlorvos	methyl paraoxon
butachlor	fenarimol
hexazinon	napropamide
bladex	prometryn
propazine	trifluralin
EPTC (s-ethyl)	metolachlor
terbutryne	bromacil
Cycloate	molinate
triadimefon	vernolate
ethoprophos (prophos)	chlorpyrifos
tetrachlorvinphos (Stirophos)	simetryn
diphenamid	fluridone
norflurazon	tricyclazole

NY/T 761-2008 Pyrethroids Pesticide Mixture 551

DRE-A50000551AC 1ml
NY/T 761-2008 Pyrethroids Pesticide Mixture 551 100 µg/mL in Acetone(‡)

dimethoate	dichlorvos
phosalone	omethoate
chlorpyrifos	methidathion
triazophos	acephate
methyl parathion	fenitrothion
malathion	methamidophos
isocarbophos	isofenphos-methyl
quinalphos	profenofos

Montana Pesticide Mixture 270

DRE-A50000270AL 1ml
Montana Pesticide Mixture 270 10 µg/mL in Acetonitrile(‡)
DRE-S50000270AL 5x1ml
Montana Pesticide Mixture 270 10 µg/mL in Acetonitrile(‡)

Abamectin	Acequinocyl
Chlormequat chloride	Cyfluthrin
Fenoxycarb	Imazalil
Paclobutrazol	Pyrethrins
Trifloxystrobin	Bifenthrin
Bifenazate	Ettoxazole
Daminozide	Myclobutanil
Imidacloprid	Spirotetramat
Spinosad	

Nevada Pesticide Mixture 694 Version 2

DRE-GA090000694AL 1ml
Nevada Pesticide Mixture Version 2 100 µg/mL in Acetonitrile(‡)(*)
DRE-GS090000694AL 5x1ml
Nevada Pesticide Mixture Version 2 100 µg/mL in Acetonitrile(‡)(*)

abamectin	acequinocyl
bifenthrin	cypermethrin (mix of isomers)
etoxazole	fenhexamid
imidacloprid	paclobutrazol (mix of isomers)
pyrethrin (mixture of isomers)	spinetoram (mixture of isomers)
Systhane TM	thiamethoxam
baythroid (mixture of isomers)	bifenazate
daminozide	dimethomorph
flonicamid	fludioxonil
pentachloronitrobenzene	piperonyl butoxide
spinosad (Spinosyn A & D)	spirotetramat
trifloxystrobin	

NY/T 761-2008 Organophosphorus Pesticide Mixture 678

DRE-A50000678AC 1ml
NY/T 761-2008 Organophosphorus Pesticide Mixture 678 1000 µg/mL in Acetone(‡)

methamidophos	parathion
isocarbophos	malathion
phosalone	terbufos
tetraethyl dithiopyrophosphate	methidathion
isofenphos-methyl	

NY/T 761-2008 Pesticide Mixture 674

DRE-A50000674AC 1ml
NY/T 761-2008 Pesticide Mixture 674 1000 µg/mL in Acetone(‡)

phorate	phorat-sulfone
trichlorfon	fenitrothion
profenofos	coumaphos
ethoprophos (prophos)	diazinon
omethoate	

Ohio Pesticide Mixture 335

DRE-A50000335AL 1ml
Ohio Pesticide Mixture 335 10-100 µg/mL in Acetonitrile(‡)
DRE-S50000335AL 5x1ml
Ohio Pesticide Mixture 335 10-100 µg/mL in Acetonitrile(‡)

Abamectin [10 µg/mL]	Aldicarb [10 µg/mL]
Daminozide [10 µg/mL]	Diazinon [100 µg/mL]
Ettoxazole [10 µg/mL]	Flonicamid [30 µg/mL]
Myclobutanil [10 µg/mL]	Paclobutrazol [10 µg/mL]
Spinosad [10 µg/mL]	Spirotetramat [100 µg/mL]
Bifenazate [20 µg/mL]	Cyfluthrin [10 µg/mL]
Dichlorvos [10 µg/mL]	Dimethoate [10 µg/mL]
Fludioxonil [10 µg/mL]	Imidacloprid [10 µg/mL]
Piperonyl butoxide [100 µg/mL]	Pyrethrins [50 µg/mL]
Thiamethoxam [20 µg/mL]	Trifloxystrobin [20 µg/mL]

Nevada Pesticide Mixture 694 Version 2

DRE-GA09000694AL 1ml
 Nevada Pesticide Mixture Version 2 100 µg/mL in Acetonitrile(±)(*
 DRE-GS09000694AL 5x1ml
 Nevada Pesticide Mixture Version 2 100 µg/mL in Acetonitrile(±)(*
 abamectin acequinocyl
 bifenthrin cypermethrin (mix of isomers)
 etoxazole fenhexamid
 imidacloprid paclobutrazol (mix of isomers)
 pyrethrin (mixture of isomers) spinetoram (mixture of isomers)
 Systhane TM thiamethoxam
 baythroid (mixture of isomers) bifenazate
 daminozide dimethomorph
 flonicamid fludioxonil
 pentachloronitrobenzene piperonyl butoxide
 spinosad (Spinosyn A & D) spirotetramat
 trifloxystrobin

NY/T 761-2008 Organophosphorus Pesticide Mixture 678

DRE-A50000678AC 1ml
 NY/T 761-2008 Organophosphorus Pesticide Mixture 678
 1000 µg/mL in Acetone(±)

methamidophos parathion
 isocarbophos malathion
 phosalone terbufos
 tetraethyl dithiopyrophosphate methidathion
 isofenphos-methyl

NY/T 761-2008 Pesticide Mixture 674

DRE-A50000674AC 1ml
 NY/T 761-2008 Pesticide Mixture 674 1000 µg/mL in Acetone(±)

phorate phorat-sulfone
 trichlorfon fenitrothion
 profenofos coumaphos
 ethoprophos (prophos) diazinon
 omethoate

Ohio Pesticide Mixture 335

DRE-A50000335AL 1ml
 Ohio Pesticide Mixture 335 10-100 µg/mL in Acetonitrile(±)
 DRE-S50000335AL 5x1ml
 Ohio Pesticide Mixture 335 10-100 µg/mL in Acetonitrile(±)

Abamectin [10 µg/mL] Aldicarb [10 µg/mL]
 Daminozide [10 µg/mL] Diazinon [100 µg/mL]
 Etoxazole [10 µg/mL] Flonicamid [30 µg/mL]
 Myclobutanil [10 µg/mL] Paclobutrazol [10 µg/mL]
 Spinosad [10 µg/mL] Spirotetramat [100 µg/mL]
 Bifenazate [20 µg/mL] Cyfluthrin [10 µg/mL]
 Dichlorvos [10 µg/mL] Dimethoate [10 µg/mL]
 Fludioxonil [10 µg/mL] Imidacloprid [10 µg/mL]
 Piperonyl butoxide [100 µg/mL] Pyrethrins [50 µg/mL]
 Thiamethoxam [20 µg/mL] Trifloxystrobin [20 µg/mL]

Oklahoma Pesticide Mixture 341

DRE-A50000341AL 1ml
 Oklahoma Pesticide Mixture 341 10 µg/mL in Acetonitrile(±)(*
 Avermectin B1 Azoxystrobin
 Bifenazate Etoxazole
 Tebuconazole Enilconazole
 Imidacloprid Malathion
 Myclobutanil Permethrin
 Spinosad Spiromesifen
 Spirotetramat

Nitrogen/Phosphorus Pesticide Mixture 925

DRE-GA09000925AC 1ml
 Nitrogen/Phosphorus Pesticide Mixture 925 100 µg/mL in Acetone(±)

alachlor atrazine
 propyzamide (pronamide) propachlor
 ametryne atraton
 prometon tebuthiuron
 butylate chlorpropham
 pebulate terbacil
 MGK-264 (mixture of isomers) phosdrinTM (mevinphos)
 dichlorvos methyl paraoxon
 butachlor fenarimol
 hexazinon napropamide
 bladex prometryn
 propazine trifluralin
 EPTC (s-ethyl) metolachlor
 terbutryne bromacil
 Cycloate molinate
 triadimefon vernolate
 ethoprophos (prophos) chlorpyrifos
 tetrachlorvinphos (Stirophos) simetryn
 diphenamid fluridone
 norflurazon tricyclazole

NY/T 761-2008 Pesticide Mixture 615

DRE-A50000615HE 1ml
 NY/T 761-2008 Pesticide Mixture 615 100 µg/mL in Hexane(±)

cypermethrin (mix of isomers) deltamethrin
 fenvalerate (mixt. of diastereoisomers) fenpropathrin
 lambda cyhalothrin baythroid (mixture of four isomers)
 bifenthrin

NY/T 761-2008 Pyrethroids Pesticide Mixture 551

DRE-A50000551AC 1ml
 NY/T 761-2008 Pyrethroids Pesticide Mixture 551 100 µg/mL in Acetone(±)

dimethoate dichlorvos
 phosalone omethoate
 chlorpyrifos methidathion
 triazophos acephate
 methyl parathion fenitrothion
 malathion methamidophos
 isocarbophos isofenphos-methyl
 quinalphos profenofos

Ohio Residual Pesticide Mixture

DRE-S50000055AL 5x1ml
 Ohio Residual Pesticide Mixture 100 µg/mL in Acetonitrile(±)(*
 Daminozide Imidacloprid
 Spinosad (mixt. of Spinosyn A and D) Flonicamid
 Pyrethrin (mixt. of isomers) Thiamethoxam
 Fludioxonil Systhane Tm
 Bifenazate Etoxazole
 Dichlorvos Aldicarb
 Dimethoate Diazinon
 Abamectin Paclobutrazol (mixt. of Stereo Isomers)
 Trifloxystrobin Piperonyl Butoxide
 Spirotetramat Baythroid (mixt. of four Isomers)

Oregon Pesticide Mixture 1

DRE-GA09000231AL 1ml
 Oregon Pesticide Mixture 1 600 µg/mL in Acetonitrile(±)
 DRE-GS09000231AL 5x1ml
 Oregon Pesticide Mixture 1 600 µg/mL in Acetonitrile(±)

abamectin spinosad (Spinosyn A & D)

Oregon Pesticide Mixture 1-100

DRE-GA09000058AL 1ml
Oregon Pesticide Mixture 1-100 100 µg/mL in Acetonitrile(±)(*)
DRE-GS09000058AL 5x1ml
Oregon Pesticide Mixture 1-100 100 µg/mL in Acetonitrile(±)(*)

abamectin	acephate
aldicarb	azoxystrobin
boscalid	carbaryl
chlorfenapyr	chlorpyrifos
cypermethrin (mix of isomers)	daminozide
acequinocyl	acetamiprid
bifenazate	bifenthrin
carbofuran	chlorantraniliprole
clofentezine	baythroid (mixture of isomers)
dichlorvos	diazinon

Oregon Pesticide Mixture 2 100x AL

DRE-GA09000256AL 1ml
Oregon Pesticide Mixture 2 100x Action Limit
100-200 µg/mL in Acetonitrile(±)(*)
DRE-GS09000256AL 5x1ml
Oregon Pesticide Mixture 2 100x Action Limit
100-200 µg/mL in Acetonitrile(±)(*)

acephate [40 µg/mL]	aldicarb [40 µg/mL]
fenpyroximate [40 µg/mL]	kresoxim methyl [40 µg/mL]
dibrom [50 µg/mL]	propiconazol (mixt. isomers) [40 µg/mL]
paclobutrazol (mixt. isomers) [40 µg/mL]	fenpropramid [40 µg/mL]
boscalid [40 µg/mL]	ethofenprox [40 µg/mL]
imidacloprid [40 µg/mL]	methomyl [40 µg/mL]
spiroxamine (mixt. isomers) [40 µg/mL]	tebuconazol (Folicur) [40 µg/mL]
abamectin [50 µg/mL]	fludioxonil [40 µg/mL]

Oregon Pesticide Mixture 2-100

DRE-GA09000059AL 1ml
Oregon Pesticide Mixture 2-100 100 µg/mL in Acetonitrile(±)
DRE-GS09000059AL 5x1ml
Oregon Pesticide Mixture 2-100 100 µg/mL in Acetonitrile(±)

dimethoate	ethoprophos (prophos)
fenpyroximate	fenpropramid
hexythiazox	imazalil
metalaxyl	methiocarb
MGK-264 - isomer b	Systhane TM
ethofenprox	fenoxycarb
flonicamid	fludioxonil
imidacloprid	kresoxim methyl
methomyl	methyl parathion
malathion	etoxazole

Oregon Pesticide Mixture 3 100x AL

DRE-GA09000257AL 1ml
Oregon Pesticide Mixture 3 100x Action Limit
100-200 µg/mL in Acetonitrile(±)(*)
DRE-GS09000257AL 5x1ml
Oregon Pesticide Mixture 3 100x Action Limit
100-200 µg/mL in Acetonitrile(±)(*)

acetamiprid	azoxystrobin
chlorpyrifos	diazinon
etoxazole	fenoxycarb
metalaxyl	methiocarb
Systhane TM	permethrin (mixture of isomers)
pyridaben	spinosad (Spinopsyn A & D)
thiamethoxam	trifloxystrobin
carbaryl	chlorantraniliprole
bifenthrin	carbofuran
dimethoate	ethoprophos (prophos)
imazalil	malathion
methyl parathion	MGK-264 - isomer b
phosmet	propoxur
spiromesifen	thiacloprid
spirotetramat	bifenazate
clofentezine	prallethrin

Oregon Pesticide Mixture 10x AL

DRE-GA09000244AL 1ml
Oregon Pesticide Mixture 10x Action Limit 2-20 µg/mL in Acetonitrile(±)(*)

abamectin [5 µg/mL]	acequinocyl [20 µg/mL]
daminozide [10 µg/mL]	dichlorvos [10 µg/mL]
flonicamid [10 µg/mL]	fludioxonil [4 µg/mL]
kresoxim methyl [4 µg/mL]	methomyl [4 µg/mL]
paclobutrazol (mixt. isomers) [4 µg/mL]	piperonyl butoxide [20 µg/mL]
tebuconazol (Folicur) [4 µg/mL]	azoxystrobin [2 µg/mL]
permethrin (mix of isomers) [2 µg/mL]	phosmet [2 µg/mL]
pyridaben [2 µg/mL]	trifloxystrobin [2 µg/mL]
diazinon [2 µg/mL]	baythroid (mixt. isomers) [10 µg/mL]
malathion [2 µg/mL]	methyl parathion [2 µg/mL]
spiromesifen [2 µg/mL]	spirotetramat [2 µg/mL]
acetamiprid [2 µg/mL]	bifenazate [2 µg/mL]
carbofuran [2 µg/mL]	chlorantraniliprole [2 µg/mL]
metalaxyl [2 µg/mL]	methiocarb [2 µg/mL]
fenoxycarb [2 µg/mL]	fenpyroximate [4 µg/mL]
aldicarb [4 µg/mL]	chlorfenapyr [10 µg/mL]
ethofenprox [4 µg/mL]	fenpropramid [4 µg/mL]
hexythiazox [10 µg/mL]	imidacloprid [4 µg/mL]
MGK-264 - isomer b [2 µg/mL]	oxamyl [10 µg/mL]
pyrethrin (mix of isomers) [10 µg/mL]	spiroxamine (mixture isomers) [4 µg/mL]
bifenthrin [2 µg/mL]	ethoprophos (prophos) [2 µg/mL]
prallethrin [2 µg/mL]	propiconazol (mixt. isomers) [4 µg/mL]
acephate [4 µg/mL]	chlorpyrifos [2 µg/mL]
cypermethrin (mixt. isomers) [10 µg/mL]	dimethoate [2 µg/mL]
Systhane TM [2 µg/mL]	spinosad (Spinopsyn A&D) [2 µg/mL]
thiacloprid [2 µg/mL]	thiamethoxam [2 µg/mL]
boscalid [4 µg/mL]	carbaryl [2 µg/mL]
clofentezine [2 µg/mL]	imazalil [2 µg/mL]
dibrom [5 µg/mL]	etoxazole [2 µg/mL]
propoxur [2 µg/mL]	

Oregon Pesticide Mixture 3

DRE-GA09000233AL 1ml
Oregon Pesticide Mixture 3 600 µg/mL in Acetonitrile(±)
DRE-GS09000233AL 5x1ml
Oregon Pesticide Mixture 3 600 µg/mL in Acetonitrile(±)

aldicarb	fenpropramid
flonicamid	hexythiazox
methiocarb	methomyl
oxamyl	pyridaben
thiacloprid	thiamethoxam

Oregon Pesticide Mixture 3-100

DRE-GA09000473AL 1ml
Oregon Pesticide Mixture 3-100 100 µg/mL in Acetonitrile(±)(*)

dibrom	oxamyl
phosmet	piperonyl butoxide
propoxur	pyrethrin (mix of isomers)
spiromesifen	spirotetramat
thiacloprid	thiamethoxam
paclobutrazol (mix of isomers)	permethrin (mix of isomers)
prallethrin	Propiconazol (mix of isomers)
pyridaben	spinosad (mix of Spinopsyn A&D)
spiroxamine	tebuconazol (Folicur)
trifloxystrobin	

Oregon Pesticide Mixture 4

DRE-GA09000234AL 1ml
Oregon Pesticide Mixture 4 600 µg/mL in Acetonitrile(±)(*)
DRE-GS09000234AL 5x1ml
Oregon Pesticide Mixture 4 600 µg/mL in Acetonitrile(±)(*)

carbaryl	carbofuran
chlorantraniliprole	clofentezine
daminozide	fenoxycarb
Imazalil	Systhane TM
paclobutrazol (mixture of stereo isomers)	Propiconazol (mixture of isomers)

Oregon Pesticide Mixture 5

DRE-GA09000235AL 1ml
Oregon Pesticide Mixture 5 600 µg/mL in Acetonitrile(±)
DRE-GS09000235AL 5x1ml
Oregon Pesticide Mixture 5 600 µg/mL in Acetonitrile(±)

bifenthrin	baythroid (mixture of isomers)
cypermethrin (mix of isomers)	ethofenprox
permethrin (mixture of isomers)	prallethrin
pyrethrin (mixture of isomers)	

Organochlorine Pesticide Mix 1

DRE-XA06080100TO 1ml
Organochlorine Pesticide Mix 1 100-600 µg/mL in Toluene

4,4'-DDD [600 µg/mL]	4,4'-DDE [200 µg/mL]
4,4'-DDT [600 µg/mL]	Aldrin [100 µg/mL]
alpha-Endosulfan [200 µg/mL]	alpha-HCH [100 µg/mL]
beta-HCH [100 µg/mL]	delta-HCH [100 µg/mL]
Dieldrin [200 µg/mL]	Endosulfan-sulfate [600 µg/mL]
Endrin [200 µg/mL]	Endrin-aldehyde [600 µg/mL]
gamma-HCH [100 µg/mL]	Heptachlor [100 µg/mL]
Heptachlor-endo-epoxide (trans-isom. A) [100 µg/mL]	

Organochlorine Pesticide Mix 3

DRE-LA06170100TH 1ml
Organochlorine Pesticide Mix 3 20 µg/mL in Toluene/Hexane
DRE-YA06170100TH 1ml
Organochlorine Pesticide Mix 3 1000 µg/mL in Toluene/Hexane

4,4'-DDD	4,4'-DDE
alpha-Endosulfan	alpha-HCH
delta-HCH	Dieldrin
Endrin-aldehyde	gamma-HCH
Methoxychlor	Aldrin
4,4'-DDT	beta-HCH
beta-Endosulfan	Endrin
Endosulfan-sulfate	Heptachlor-endo-epoxide (trans-isom. A)
Heptachlor	

Organochlorine Pesticide Mixture 372

DRE-GA09000372TH 1ml
Organochlorine Pesticide Mixture 372 200 µg/mL in Toluene:Hexane(±)

aldrin	α-BHC (alpha-HCH)
δ-BHC (delta-HCH)	cis-chlordane
p,p'-DDE	p,p'-DDT
endosulfan II	endosulfan sulfate
endrin ketone	heptachlor
β-BHC (beta-HCH)	γ-BHC (gamma-HCH, Lindane)
trans-chlordane	p,p'-DDD
dieldrin	endosulfan I
endrin	endrin aldehyde
heptachlor epoxide isomer B	methoxychlor

Oregon Pesticide Mixture 476

DRE-A50000476AL 1ml
Oregon Pesticide Mixture 476 600 µg/mL in Acetonitrile(±)
DRE-S50000476AL 5x1ml
Oregon Pesticide Mixture 476 600 µg/mL in Acetonitrile(±)

(E)-Fenpyroximate	Acequinocyl
Bifenazate	Boscalid
Fludioxonil	Imidacloprid
MGK 264	Piperonyl butoxide
Spiroxamine	Trifloxystrobin
Acetamiprid	Azoxystrobin
Chlorfenapyr	Etiozazole
Kresoxim-methyl	Metaxalyl
Spiromesifen	Spirotetramat

Oregon Pesticide Mixture 662 100x AL

DRE-A50000662AL 1ml
Oregon Pesticide Mixture 662 100x Action Limit 100-200 µg/mL in Acetonitrile(±)

Acequinocyl [200 µg/mL]	Chlorfenapyr [100 µg/mL]
Cyfluthrin [100 µg/mL]	Cypermethrin (technical) [100 µg/mL]
Daminozide [100 µg/mL]	Dichlorvos [100 µg/mL]
Fonicamid [100 µg/mL]	Hexythiazox [100 µg/mL]
Oxamyl [100 µg/mL]	Piperonyl butoxide [200 µg/mL]
Pyrethrins [100 µg/mL]	

Organochlorine Pesticide Mix 2

DRE-YA08080100TH 1ml
Organochlorine Pesticide Mix 2 2000 µg/mL in Toluene/Hexane(±)

4,4'-DDD	4,4'-DDE
alpha-Endosulfan	alpha-HCH
delta-HCH	Dieldrin
Endrin-aldehyde	Endrin-ketone
Heptachlor-endo-epoxide (trans-isom. A) Methoxychlor	Aldrin
4,4'-DDT	beta-HCH
beta-Endosulfan	Endrin
Endosulfan-sulfate	Heptachlor
gamma-HCH	

Organochlorine Pesticide Mixture 57

DRE-GS09000057DI 5x1ml
Organochlorine Pesticide Mixture 57 200 µg/mL in Dichloromethane(±)

a-BHC	b-BHC
aldrin	heptachlor
trans-chlordane	p,p'-DDD
dieldrin	endrin
endosulfan I	endosulfan II
d-BHC	g-BHC
heptachlor epoxide isomer B	cis-chlordane
p,p'-DDE	p,p'-DDT
endrin aldehyde	endrin ketone
endosulfan sulfate	methoxychlor

Organochlorine Pesticide Mixture 373

DRE-GA09000373TH 1ml
Organochlorine Pesticide Mixture 373 1000 µg/mL in Toluene:Hexane(±)(*)

aldrin	α-BHC (alpha-HCH)
δ-BHC (delta-HCH)	cis-chlordane
p,p'-DDE	p,p'-DDT
endosulfan II	endosulfan sulfate
endrin ketone	heptachlor
β-BHC (beta-HCH)	γ-BHC (gamma-HCH, Lindane)
trans-chlordane	p,p'-DDD
dieldrin	endosulfan I
endrin	endrin aldehyde
heptachlor epoxide isomer B	methoxychlor

Organochlorine Pesticides Mixture 921

DRE-GA09000921AC 1ml
Organochlorine Pesticides Mixture 921 2000 µg/mL in Acetone(‡)

endosulfan I	endosulfan II
p,p'-DDE	p,p'-DDT
endrin aldehyde	endrin ketone
b-BHC	d-BHC
heptachlor	heptachlor epoxide isomer B
endosulfan sulfate	p,p'-DDD
dieldrin	endrin
methoxychlor	a-BHC
g-BHC	aldrin

Organochlorine Pesticide Mixture 946

DRE-GA09000946AC 1ml
Organochlorine Pesticide Mixture 946 100-200 µg/mL in Acetone(‡)

p,p'-DDT [100 µg/mL]	p,p'-DDD [100 µg/mL]
g-BHC [100 µg/mL]	dieldrin [100 µg/mL]
endosulfan sulfate [100 µg/mL]	endrin [100 µg/mL]
b-BHC [100 µg/mL]	d-BHC [100 µg/mL]
cis-chlordane [100 µg/mL]	trans-chlordane [100 µg/mL]
etridiazole [100 µg/mL]	alachlor [100 µg/mL]
simazine [100 µg/mL]	permethrin (mix of isomers) [200 µg/mL]
chloroneb [100 µg/mL]	a-BHC [100 µg/mL]
p,p'-DDE [100 µg/mL]	endosulfan II [100 µg/mL]
methoxychlor [100 µg/mL]	aldrin [100 µg/mL]
endrin aldehyde [100 µg/mL]	heptachlor epoxide isom B [100µg/mL]
heptachlor [100 µg/mL]	trans-nonachlor [100 µg/mL]
endosulfan I [100 µg/mL]	chlorobenzilate [100 µg/mL]
atrazine [100 µg/mL]	chlorothalonil [100 µg/mL]
chlorthal-dimethyl (dacthal) [100 µg/mL]	

Organochlorine Pesticide Mixture 952

DRE-GA09000952TH 1ml
Organochlorine Pesticide Mixture 952 1000 µg/mL in Toluene:Hexane(‡)

methoxychlor	a-BHC
g-BHC	aldrin
cis-chlordane	trans-chlordane
p,p'-DDT	dieldrin
endrin ketone	endosulfan I
b-BHC	d-BHC
heptachlor	heptachlor epoxide isomer B
p,p'-DDD	p,p'-DDE
endrin	endrin aldehyde
endosulfan II	endosulfan sulfate

Organochlorine Pesticide Mixture 1012

DRE-GA09001012AC 1ml
Organochlorine Pesticide Mixture 1012 2000 µg/mL in Acetone(‡)

p,p'-DDD	p,p'-DDE
endrin	a-BHC
g-BHC	aldrin
endosulfan I	endosulfan II
endrin aldehyde	dieldrin
p,p'-DDT	d-BHC
b-BHC	heptachlor epoxide isomer B
heptachlor	endosulfan sulfate
methoxychlor	

Organochlorine Pesticides Internal Standards Mixture 135 for HJ 835-2017, HJ 900, HJ 912-2017

DRE-A50000135AH 1ml
HJ 835-2017, HJ 900, HJ 912-2017 Organochlorine Pesticides Internal Standards Mixture 135 1000 µg/mL in Acetone:n-Hexane(‡)(*)
Phenanthrene D10 Quintozene

Organochlorine Pesticides Mixture 942

DRE-GA09000942MB 1ml
Organochlorine Pesticides Mixture 942 1000 µg/mL in Methyl tert-butyl ether(‡)

p,p'-DDD	p,p'-DDE
endrin	a-BHC
g-BHC	aldrin
endosulfan I	endosulfan II
endrin aldehyde	dieldrin
p,p'-DDT	d-BHC
b-BHC	heptachlor epoxide isomer B
heptachlor	endosulfan sulfate
methoxychlor	

Organochlorine Pesticide Mixture 951

DRE-GA09000951ME 1ml
Organochlorine Pesticide Mixture 951 20-100 µg/mL in Methanol(‡)(*)

p,p'-DDD [100 µg/mL]	p,p'-DDT [100 µg/mL]
endosulfan II [100 µg/mL]	endosulfan sulfate [100 µg/mL]
d-BHC [20 µg/mL]	g-BHC [20 µg/mL]
heptachlor epoxide isomer B [20 µg/mL]	p,p'-DDE [20 µg/mL]
methoxychlor [20 µg/mL]	endrin [100 µg/mL]
dieldrin [20 µg/mL]	b-BHC [20 µg/mL]
a-BHC [20 µg/mL]	heptachlor [20 µg/mL]
aldrin [20 µg/mL]	endosulfan I [20 µg/mL]
endrin aldehyde [20 µg/mL]	

Organochlorine Pesticide Mixture 1008

DRE-GA09001008TH 1ml
Organochlorine Pesticide Mixture 1008 2000 µg/mL in Toluene:Hexane(‡)

aldrin	a-BHC
d-BHC	p,p'-DDD
dieldrin	endosulfan I
endrin	endrin aldehyde
b-BHC	heptachlor
p,p'-DDE	p,p'-DDT
endosulfan II	endosulfan sulfate
g-BHC	heptachlor epoxide isomer A

Organochlorine Pesticides Decomposition Mixture

DRE-S50000285EA 5x1ml
Organochlorine Pesticides Decomposition Mixture 100 µg/mL in Ethyl acetate(‡)
4,4'-DDT Endrin

Organochlorine Pesticides Mixture 105 for HJ 835-2017, HJ 900, HJ 901, HJ 912-2017

DRE-A50000105AH 1ml
HJ 835-2017, HJ 900, HJ 901, HJ 912-2017 Organochlorine Pesticides Mixture 105 1000 µg/mL in Acetone:n-Hexane(‡)(*)

Aldrin	beta-Endosulfan
alpha-HCH	beta-HCH
Heptachlor	2,4'-DDT
4,4'-DDD	Methoxychlor
Dieldrin	Endrin-aldehyde
Endosulfan-sulfate	Mirex
alpha-Endosulfan	Hexachlorobenzene
delta-HCH	gamma-HCH
4,4'-DDT	4,4'-DDE
Dicofol	Endrin
Endrin-ketone	Heptachlor-exo-epoxide (cis-, isomer B)
Cis-Chlordane (Alpha Isomer)	Trans-Chlordane (Gamma Isomer)

Organochlorine Pesticides Mixture 110 for GB 2763

DRE-A50000111TH 1ml
 GB 2763 Organochlorine Pesticides Mixture 110
 20-100 µg/mL in Toluene:n-Hexane(‡)

Aldrin [100 µg/mL]	beta-Endosulfan [100 µg/mL]
beta-HCH [100 µg/mL]	delta-HCH [100 µg/mL]
2,4'-DDT [100 µg/mL]	4,4'-DDT [100 µg/mL]
Dieldrin [100 µg/mL]	Heptachlor-exo-epoxide (B) [50 µg/mL]
Heptachlor-endo-epoxide (A) [50 µg/mL]	cis-Chlordane (alpha Isomer) [20 µg/mL]
alpha-Endosulfan [100 µg/mL]	alpha-HCH [100 µg/mL]
gamma-HCH [100 µg/mL]	Heptachlor [50 µg/mL]
4,4'-DDE [100 µg/mL]	4,4'-DDD [100 µg/mL]
Endosulfan-sulfate [100 µg/mL]	oxy-Chlordane [20 µg/mL]
trans-Chlordane (gamma) [20 µg/mL]	

Organochlorine Pesticides Mixture 302 for HJ 835-2017, HJ 835-2017, HJ 900-2017, HJ 901-2017, HJ 912-2017, HJ 921-2017

DRE-A50000302AH 1ml
 HJ 835-2017, HJ 835-2017, HJ 900-2017, HJ 901-2017,
 HJ 912-2017, HJ 921-2017 Organochlorine
 Pesticides Mixture 302 1000 µg/mL in n-Hexane/Acetone(‡)

Aldrin	Cis-Chlordane (Alpha Isomer)
4,4'-DDE	2,4'-DDT
Endosulfan-sulfate	alpha-Endosulfan
Endrin-aldehyde	Endrin-ketone
delta-HCH	gamma-HCH
Hexachlorobenzene	Methoxychlor
trans-Chlordane (Gamma Isomer)	4,4'-DDD
4,4'-DDT	Dieldrin
beta-Endosulfan	Endrin
alpha-HCH	beta-HCH
Heptachlor	Heptachlor-exo-epoxide (isom. B)
Mirex	

Organochlorine Pesticides Substitutes Mixture 129 for HJ 835-2017, HJ 912-2017

DRE-A50000129AH 1ml
 HJ 835-2017, HJ 912-2017 Organochlorine Pesticides Substitutes Mixture 129
 1000 µg/mL in Acetone:n-Hexane(‡)

2,4,5,6-Tetrachloro-m-xylene Dibutyl chlorendate

Organohalide Pesticide Mix 2

DRE-XA05050200ME 1ml
 Organohalide Pesticide Mix 2 100 µg/mL in Methanol(‡)

Alachlor	Aldrin
Dieldrin	Endrin
gamma-HCH	Heptachlor
Heptachlor-endo-epoxide (trans-isom. A)	Heptachlor-exo-epoxide (cis-isomer B)
Hexachlorobenzene	Methoxychlor

Organophosphorous Pesticide Mix 1

DRE-XA06140100AC 1ml
 Organophosphorous Pesticide Mix 1 200 µg/mL in Acetone

Azinphos-methyl	Demeton (O+S)
Diazinon	Disulfoton
Ethion	Malathion
Parathion-ethyl	Parathion-methyl

Organochlorine Pesticides Mixture 122 for GB/T 5750.9-2006, GB/T 14848-2017

DRE-A50000122TO 1ml
 GB/T 5750.9-2006, GB/T 14848-2017 Organochlorine Pesticides Mixture 122
 100 µg/mL in Toluene(‡)

alpha-HCH	beta-HCH
delta-HCH	gamma-HCH
2,4'-DDT	4,4'-DDT
4,4'-DDE	4,4'-DDD

Organochlorine Pesticides Mixture 637

DRE-A50000637HE 1ml
 Organochlorine Pesticides Mixture 637 500 µg/mL in Hexane(‡)

endosulfan I	endosulfan II
d-BHC	g-BHC
heptachlor epoxide isomer B	isodrin
p,p'-DDD	p,p'-DDE
hexachlorobenzene	o,p'-DDD
mirex	b-BHC
a-BHC	heptachlor
aldrin	p,p'-DDT
methoxychlor	endrin
dieldrin	o,p'-DDT
o,p'-DDE	

Organohalide Pesticide Mix 1

DRE-XA05050100AC 1ml
 Organohalide Pesticide Mix 1 100 µg/mL in Acetone

Alachlor	Aldrin
Dieldrin	Endrin
Heptachlor-endo-epoxide (trans-isom. A)	Heptachlor-exo-epoxide (cis-isomer B)
Methoxychlor	Simazine
Atrazine	cis-Nonachlor
gamma-HCH	Heptachlor
Hexachlorobenzene	Hexachlorocyclopentadiene
trans-Chlordane (gamma)	trans-Nonachlor

Organometallic Butyltin Chloride Mixture

DRE-A50000280DI 1ml
 Organometallic Butyltin Chloride Mixture 2000 µg/mL in Dichloromethane(‡)

Butyltin trichloride	Dibutyltin dichloride
Tetrabutyltin	Tributyltin chloride

Organophosphorous Pesticides Mix 4

DRE-XA06220100CY 1ml
 Organophosphorous Pesticides Mix 4 100 µg/mL in Cyclohexane

Azinphos-methyl	Coumaphos
Demeton (O+S)	Disulfoton
Fensulfothion	Fenthion
Phorate	Prothiophos
Sulprofos	Trichloronate

Organophosphorus Pesticide Mixture 682

DRE-A50000682CY 1ml
Organophosphorus Pesticide Mixture 682 100 µg/mL in Cyclohexane(‡)

dichlorvos	disulfoton
ethion	fenitrothion
malathion	methidathion
phosdrinTM (mevinphos)	parathion
methyl parathion	phorate
phosphamidon	triazophos
monocrotophos	

Organophosphorus Pesticides Mixture 573

DRE-A50000573ME 1ml
Organophosphorus Pesticides Mixture 573 100 µg/mL in Methanol(‡)

omethoate	isocarbophos
isofenphos-methyl	chlorpyrifos
methyl parathion	parathion
monocrotophos	ethoprophos (prophos)
acephate	isazophos
triazophos	dimethoate
phoxim	methidathion
phorate	dichlorvos
methamidophos	

Organophosphorus Pesticides Mixture 580

DRE-A50000580AH 1ml
Organophosphorus Pesticides Mixture 580 400 µg/mL in Acetone:Hexane(‡)(*)

trichloronate	sulprofos
fensulfotthion	tetrachlorvinphos (Rabon)
Guthion®	dichlorvos
malathion	tetraethyl dithiopyrophosphate
chlorpyrifos	tributylphosphoro-trithioite (Merphos)
thionazine (zinophos)	phorate
monocrotophos	phosdrinTM (mevinphos)
propazine	Tokuthion®
anilazine	fenchlorphos
coumaphos	Demeton O&S
dibrom	diazinon
simazine	tetraethyl pyrophosphate
methyl parathion	EPN
famphur	O,O,O-triethylphosphorothioate
disulfoton	parathion
ethoprophos (prophos)	atrazine
dimethoate	fenthion

Organophosphorus Pesticides Mixture 597

DRE-A50000597EA 1ml
Organophosphorus Pesticides Mixture 597 100 µg/mL in Ethyl acetate(‡)

methyl parathion	dimethoate
omethoate	methamidophos
monocrotophos	diazinon
ethion	malathion
methidathion	dichlorvos
acephate	parathion

Organophosphorus Pesticides Mixture 624

DRE-A50000624MB 1ml
Organophosphorus Pesticides Mixture 624
1000 µg/mL in Methyl tert Butyl Ether(‡)

methyl parathion	parathion
malathion	dimethoate
Demeton O&S	trichlorfon
methamidophos	dichlorvos

Organophosphorus Pesticides Mixture 665

DRE-A50000665AC 1ml
Organophosphorus Pesticides Mixture 665 100 µg/mL in Acetone(‡)

methamidophos	acephate
parathion	isocarbophos
methidathion	triazophos
monocrotophos	methyl parathion

Organophosphorus Pesticides Mixture 988

DRE-GA09000988AC 1ml
Organophosphorus Pesticides Mixture 988 200 µg/mL in Acetone(‡)

Guthion®	Demeton O&S
diazinon	disulfoton
ethion	malathion
parathion	methyl parathion

Pennsylvania Pesticide Mixture

DRE-A50000333AL 1ml Pennsylvania Pesticide Mixture 333 10 µg/mL in Acetonitrile(‡)(*)
DRE-A50000334AL 1ml Pennsylvania Pesticide Mixture 334 100 µg/mL in Acetonitrile(‡)(*)

Abamectin	Acephate	Malathion	Metalaxyl	Imidacloprid	Kresoxim-methyl
Aldicarb	Azoxystrobin	MGK 264	Myclobutanil	Methiocarb	Methomyl
Boscalid	Captan	Acequinocyl	Acetamiprid	Naled	Oxamyl
Chlorantraniliprole	Chlorfenapyr	Bifenazate	Bifenthrin	Pacllobutrazol	Parathion-methyl
Cyfluthrin	Cypermethrin (technical)	Carbaryl	Carbofuran	Piperonyl butoxide	Prallethrin
Dichlorvos	Dimethoate	Chlorpyrifos	Clofentezine	Pyridaben	Spinetoram
Etofenprox	Etoxazole	Daminozide	Diazinon	Spirotetramat	Spiroxamine
Fenpyroximate (E/Z Mix)	Fipronil	Dimethomorph	Ethoprophos	Thiamethoxam	Trifloxystrobin
Hexythiazox	Imazalil	Fenhexamid	Fenoxycarb	Permethrin	Propiconazole
Spinosad	Tebuconazole	Fonicamid	Fludioxonil	Phosmet	Propoxur
Spiromesifen	Thiacloprid				

Pesticide-Mix 2

DRE-L18000002CY 10ml
Pesticide-Mix 2 10 µg/mL in Cyclohexane

Diazinon	Ethion
Malathion	Parathion-ethyl
Parathion-methyl	

Pesticide Mix 5

DRE-L18000005CY 10ml
Pesticide Mix 5 10 µg/mL in Cyclohexane(‡)

2,4'-DDD	2,4'-DDE
2,4'-DDT	4,4'-DDD
4,4'-DDE	4,4'-DDT
Dieldrin	Endrin

Pesticide Mix 7

DRE-L1800007CY 10ml
 Pesticide Mix 7 10 µg/mL in Cyclohexane(±)
 alpha-HCH beta-HCH
 gamma-HCH delta-HCH

Pesticide-Mix 10

DRE-Z18000010IO 10ml
 Pesticide-Mix 10 1 µg/mL in Isooctane
 2,4'-DDE 4,4'-DDD
 2,4'-DDE 4,4'-DDE
 Dieldrin Endrin
 alpha-HCH beta-HCH
 gamma-HCH delta-HCH
 Heptachlor-endo-epoxide (trans-isom. A) Hexachlorobenzene
 Methoxychlor Mirex

Pesticide Mix 13

DRE-LA18000013CY 1ml
 Pesticide-Mix 13 10 µg/mL in Cyclohexane(±)
 DRE-L18000013CY 10ml
 Pesticide Mix 13 10 µg/mL in Cyclohexane(±)
 DRE-LA18000013TO 1ml
 Pesticide-Mix 13 10 µg/mL in Toluene(±)
 DRE-L18000013TO 10ml
 Pesticide Mix 13 10 µg/mL in Toluene(±)

2,4'-DDD 2,4'-DDE
 4,4'-DDE 4,4'-DDT
 alpha-HCH beta-Endosulfan
 delta-HCH Dieldrin
 gamma-HCH Heptachlor
 Hexachlorobenzene Isodrin
 oxy-Chlordane PCB 101
 PCB 180 PCB 28
 2,4'-DDT 4,4'-DDD
 Aldrin alpha-Endosulfan
 beta-HCH cis-Chlordane (alpha)
 Endrin epsilon-HCH
 Heptachlor-endo-epoxide (trans-isom. A) Heptachlor-exo-epoxide (cis-isomer B)
 Methoxychlor Mirex
 PCB 138 PCB 153
 PCB 52 trans-Chlordane (gamma)

Pesticide-Mix 17

DRE-L18000017CY 10ml
 Pesticide-Mix 17 10 µg/mL in Cyclohexane(±)

1,2,4,5-Tetrachlorobenzene 1,2,4-Trichlorobenzene
 4,4'-DDE 4,4'-DDT
 alpha-HCH beta-Endosulfan
 Endrin gamma-HCH
 Hexachlorobenzene Methoxychlor
 2,4'-DDT 4,4'-DDD
 Aldrin alpha-Endosulfan
 beta-HCH Dieldrin
 Heptachlor Heptachlor-endo-epoxide (trans-isom. A)
 Pentachlorobenzene Quintozene

Pesticide-Mix 20

DRE-LA18000020CY 1ml
 Pesticide-Mix 20 10 µg/mL in Cyclohexane(±)
 DRE-L18000020CY 10ml
 Pesticide-Mix 20 10 µg/mL in Cyclohexane(±)

2,4'-DDD 2,4'-DDE
 4,4'-DDE 4,4'-DDT
 alpha-HCH beta-Endosulfan
 Endrin gamma-HCH
 Hexachlorobenzene Isobenzan
 2,4'-DDT 4,4'-DDD
 Aldrin alpha-Endosulfan
 beta-HCH Dieldrin
 Heptachlor Heptachlor-exo-epoxide (cis-isomer B)
 Isodrin

Pesticide-Mix 8

DRE-LA18000008CY 1ml
 Pesticide-Mix 8 10 µg/mL in Cyclohexane(±)
 alpha-HCH beta-HCH
 gamma-HCH delta-HCH
 epsilon-HCH

Pesticide Mix 11

DRE-L18000011CY 10ml
 Pesticide Mix 11 10 µg/mL in Cyclohexane(±)
 alpha-HCH beta-HCH
 gamma-HCH delta-HCH
 Hexachlorobenzene

Pesticide Mix 14

DRE-LA18000014CY 1ml
 Pesticide-Mix 14 10 µg/mL in Cyclohexane(±)
 DRE-L18000014CY 10ml
 Pesticide Mix 14 10 µg/mL in Cyclohexane(±)

Aldrin 4,4'-DDD
 4,4'-DDT Dieldrin
 Endrin alpha-HCH
 Heptachlor Heptachlor-endo-epoxide (trans-isom. A)
 4,4'-DDE 2,4'-DDT
 alpha-Endosulfan beta-Endosulfan
 beta-HCH gamma-HCH
 Hexachlorobenzene Methoxychlor

Pesticide-Mix 15

DRE-LA18000015AC 1ml
 Pesticide-Mix 15 10 µg/mL in Acetone(±)

Atrazine Atrazine-desethyl
 Atrazine-desisopropyl Propazine
 Simazine Terbutylazine

Pesticide-Mix 18

DRE-XA18000018AL 1ml
 Pesticide-Mix 18 100 µg/mL in Acetonitrile(±)

Atrazine Atrazine-desethyl
 Atrazine-desisopropyl Propazine
 Sebutylazine Simazine
 Terbutylazine

Pesticide-Mix 31

DRE-L18000031AC 10ml
 Pesticide-Mix 31 10 µg/mL in Acetone

Atrazine Atrazine-desethyl
 Metazachlor Metolachlor
 Simazine Terbutylazine

Pesticide-Mix 32

DRE-L18000032CY 10ml
Pesticide-Mix 32 10 µg/mL in Cyclohexane

2,4,5-T-methyl ester	2,4-D-methyl ester
Dichlorprop-methyl ester	MCPA-methyl ester
MCPB-methyl ester	Mecoprop-methyl ester

Pesticide-Mix 34

DRE-L18000034AL 10ml
Pesticide-Mix 34 10 µg/mL in Acetonitrile(‡)

Atrazine	Atrazine-desethyl
Chloroxuron	Chlorpropham
Diuron	Fenuron
Metamitron	Metazachlor
Metolachlor	Metoxuron
Prometryn	Propazine
Simazine	Terbuthylazine
Atrazine-desethyl-desisopropyl	Chlorotoluron
Crimidine	Cyanazine
Isoproturon	Linuron
Methabenzthiazuron	Metobromuron
Metribuzin	Monolinuron
Propham	Sebuthylazine
Terbutryn	

Pesticide-Mix 41

DRE-L18000041AL 10ml
Pesticide-Mix 41 10 µg/mL in Acetonitrile

Atrazine	Atrazine-desethyl
Atrazine-desisopropyl	Propazine
Sebuthylazine	Simazine
Terbuthylazine	Terbuthylazine-desethyl

Pesticide Mix 51

DRE-L18000051AL 10ml
Pesticide Mix 51 10 µg/mL in Acetonitrile

Atrazine	Atrazine-desethyl
Atrazine-desisopropyl	Cyanazine
Prometryn	Propazine
Sebuthylazine	Simazine
Terbuthylazine	Terbutryn

Pesticide-Mix 60

DRE-LA18000060AC 1ml
Pesticide-Mix 60 10 µg/mL in Acetone

Atrazine	Atrazine-desethyl
Atrazine-desisopropyl	Propazine
Sebuthylazine	Sebuthylazine-desethyl
Simazine	Terbuthylazine
Terbuthylazine-desethyl	

Pesticide-Mix 64

DRE-L18000064CY 10ml
Pesticide-Mix 64 10 µg/mL in Cyclohexane(*)

Chlorpyrifos	Chlorpyrifos-methyl
Diazinon	Dichlorvos
Fenitrothion	Malathion
Methacrifos	Phosphamidon
Pirimiphos-methyl	

Pesticide-Mix 33

DRE-LA18000033IO 1ml
Pesticide-Mix 33 10 µg/mL in Isooctane

4,4'-DDD	4,4'-DDE
alpha-Endosulfan	alpha-HCH
Dieldrin	Endrin
Heptachlor-endo-epoxide (trans-isom. A)	Heptachlor-exo-epoxide (cis-isomer B)
Isobenzan	Isodrin
PCB 138	PCB 153
PCB 52	Quintozene
4,4'-DDT	Aldrin
beta-HCH	Dichlobenil
gamma-HCH	Heptachlor
Hexachloro-1,3-butadiene	Hexachlorobenzene
PCB 101	PCB 118
PCB 180	PCB 28

Pesticide-Mix 40

DRE-LA18000040AL 1ml
Pesticide-Mix 40 10 µg/mL in Acetonitrile(‡)

DRE-L18000040AL 10ml
Pesticide-Mix 40 10 µg/mL in Acetonitrile

2,4-D	2,4-DB
Dichlorprop	MCPA
MCPB	Mecoprop

Pesticide-Mix 44

DRE-LA18000044AL 1ml
Pesticide-Mix 44 10 µg/mL in Acetonitrile

DRE-L18000044AL 10ml
Pesticide-Mix 44 10 µg/mL in Acetonitrile(‡)

Atrazine	Atrazine-desethyl
Diuron	Hexazinone
Metazachlor	Methabenzthiazuron
Metoxuron	Monolinuron
Terbuthylazine	Cyanazine
Chlorotoluron	Linuron
Isoproturon	Metolachlor
Metobromuron	Simazine
Sebuthylazine	

Pesticide-Mix 56

DRE-L18000056AL 10ml
Pesticide-Mix 56 10 µg/mL in Acetonitrile

Atrazine	Atrazine-desethyl
Atrazine-desisopropyl	Chlorotoluron
Diuron	Isoproturon
Propazine	Simazine

Pesticide-Mix 62

DRE-L18000062CY 10ml
Pesticide-Mix 62 10 µg/mL in Cyclohexane(‡)

Cypermethrin (technical)	Deltamethrin
Fenprothrin	Fenvalerate
Permethrin	Resmethrin

Pesticide-Mix 71

DRE-L18000071CY 10ml
Pesticide-Mix 71 10 µg/mL in Cyclohexane(‡)

2,4'-DDD	2,4'-DDE
4,4'-DDE	4,4'-DDT
alpha-HCH	beta-Endosulfan
Dieldrin	Endrin
Heptachlor-endo-epoxide (trans-isom. A)	Heptachlor-exo-epoxide (cis-isomer B)
Quintozene	Tecnazene
2,4'-DDT	4,4'-DDD
Aldrin	alpha-Endosulfan
beta-HCH	cis-Chlordane (alpha)
gamma-HCH	Heptachlor
Hexachlorobenzene	Methoxychlor
trans-Chlordane (gamma)	

Pesticide-Mix 95

DRE-YA18000095EA 1ml
Pesticide-Mix 95 1000 µg/mL in Ethyl acetate

2,3,6-Trichlorobenzoic acid	2,4,5-T
2,4-D	2,4-DB
Bentazone	Bromoxynil
Dicamba	Dichlorprop
loxynil	MCPA
MCPB	Mecoprop
Triclopyr	

Pesticide-Mix 102

DRE-LA18000102AL 1ml
Pesticide-Mix 102 50 µg/mL in Acetonitrile(±)

2,4'-DDT	4,4'-DDD
Aldrin	alpha-Endosulfan
beta-HCH	Chlorpyrifos
gamma-HCH	Heptachlor
Metolachlor	Alachlor
4,4'-DDT	beta-Endosulfan
alpha-HCH	Endrin
Dieldrin	Hexachlorobenzene
Heptachlor-endo-epoxide (trans-isom. A)	

Pesticide-Mix 118

DRE-L18000118CY 10ml
Pesticide-Mix 118 10 µg/mL in Cyclohexane(±)

Cyfluthrin	Cypermethrin (technical)
Deltamethrin	Fenvalerate
Permethrin	Phenothrin
Tetramethrin	

Pesticide-Mix 128

DRE-L18000128AL 10ml
Pesticide-Mix 128 20 µg/mL in Acetonitrile(±)

Clofibrac acid	2,4,5-T
Bentazone	Bromacil
Dichlorprop	Fenoprop
MCPA	MCPB
Triclopyr	2,4-DB
2,4-D	Dicamba
Clopyralid	Haloxypop (free acid)
Fluazifop (free acid)	Picloram
Mecoprop	

Pesticide-Mix 134

DRE-L18000134AL 10ml
Pesticide-Mix 134 10 µg/mL in Acetonitrile

Chlorotoluron	Diuron
Isoproturon	Linuron

Pesticide-Mix 164

DRE-L18000164CY 10ml
Pesticide-Mix 164 10 µg/mL in Cyclohexane(±)

2,4'-DDD	4,4'-DDD
2,4'-DDE	4,4'-DDE
2,4'-DDT	4,4'-DDT

Pesticide-Mix 101

DRE-LA18000101AL 1ml
Pesticide-Mix 101 50 µg/mL in Acetonitrile

Alachlor	Ametryn
Azinphos-methyl	Chlorpyrifos
Metolachlor	Molinate
Pendimethalin	Pirimicarb
Simazine	Terbutylazine
Atrazine	Azinphos-ethyl
Diazinon	Malathion
Parathion-ethyl	Parathion-methyl
Prometryn	Propazine
Terbutryn	Trifluralin

Pesticide-Mix 114

DRE-XA18000114IO 1ml
Pesticide-Mix 114 10-200 µg/mL in Isooctane

Bifenthrin [10 µg/mL]	Cyfluthrin [40 µg/mL]
Cypermethrin (technical) [40 µg/mL]	Deltamethrin [40 µg/mL]
Fenvalerate [30 µg/mL]	Permethrin [50 µg/mL]
Phenothrin [200 µg/mL]	

Pesticide-Mix 123

DRE-LA18000123AL 1ml
Pesticide-Mix 123 10-20 µg/mL in Acetonitrile

Alachlor [20 µg/mL]	Ametryn [10 µg/mL]
Atrazine [10 µg/mL]	Atrazine-desethyl [10 µg/mL]
Cyanazine [10 µg/mL]	Pendimethalin [20 µg/mL]
Prometryn [10 µg/mL]	Propazine [10 µg/mL]
Simazine [10 µg/mL]	Terbutylazine [10 µg/mL]
Terbutryn [10 µg/mL]	Trifluralin [10 µg/mL]

Pesticide-Mix 129

DRE-L18000129AL 10ml
Pesticide-Mix 129 20-40 µg/mL in Acetonitrile

Chlorotoluron [20 µg/mL]	Diuron [20 µg/mL]
Fenuron [20 µg/mL]	Fluometuron [20 µg/mL]
Isoproturon [20 µg/mL]	Linuron [20 µg/mL]
Metobromuron [20 µg/mL]	Metoxuron [40 µg/mL]
Monolinuron [20 µg/mL]	Monuron [20 µg/mL]

Pesticide-Mix 163

DRE-YA18000163TH 1ml
Pesticide-Mix 163 2000 µg/mL in Toluene/Hexane

4,4'-DDD	4,4'-DDE
alpha-Endosulfan	alpha-HCH
delta-HCH	Dieldrin
gamma-HCH	Heptachlor
4,4'-DDT	Aldrin
beta-Endosulfan	beta-HCH
Endosulfan-sulfate	Endrin
Heptachlor-endo-epoxide (trans-isom. A)	Methoxychlor

Pesticide-Mix 167

DRE-L18000167TO 10ml
Pesticide-Mix 167 10 µg/mL in Toluene(‡)

Azinphos-ethyl	Azinphos-methyl
Carbophenothion	Chlorfenvinphos
Diazinon	Dichlofenthion
Etrimfos	Fenchlorphos
Methacrifos	Methidathion
Propetamphos	Sulfotep
Bromophos-ethyl	Bromophos-methyl
Chlorpyrifos	Chlorpyrifos-methyl
Dichlorvos	Ethion
Fonofos	Malathion
Pirimiphos-ethyl	Pirimiphos-methyl
Tetrachlorvinphos	

Pesticide-Mix 192

DRE-L18000192IO 10ml
Pesticide-Mix 192 5-10 µg/mL in Isooctane

Bifenthrin [5 µg/mL]	Cyfluthrin [10 µg/mL]
Flucythrinate [10 µg/mL]	lambda-Cyhalothrin [5 µg/mL]
tau-Fluvalinate [10 µg/mL]	

Pesticide Mixture 200

DRE-GA09000200EA 1ml
Pesticide Mixture 200 100 µg/mL in Ethyl acetate(‡)

b-BHC	d-BHC
o,p'-DDE	pentachloronitrobenzene
a-BHC	hexachlorobenzene
aldrin	heptachlor epoxide isomer B
dieldrin	endrin
methoxychlor	p,p'-DDT
o,p'-DDD	o,p'-DDT
endosulfan I	heptachlor epoxide isomer A
mirex	oxychlordane
g-BHC	heptachlor
trans-chlordane	cis-chlordane
endosulfan II	endosulfan sulfate
p,p'-DDD	p,p'-DDE

Pesticide Mixture 205

DRE-GH09000205AL 10x1ml
Pesticide Mixture 205 100 µg/mL in Acetonitrile(‡)(*)

lufenuron	metalaxyl
methoxyfenozide	novaluron
phenothoate	prochloraz
pyraclostrobin	pyrimethanil
tebuconazol (Folicur)	tebufenozide
Tolfenpyrade	triazophos
malathion	methomyl
methidathion	oxamyl
omethoate	propargite
profenofos	spinosad (Spinosyn A & D)
spinetoram (mixture of isomers)	thiacloprid
thiabendazole	imidacloprid
trifloxystrobin	

Pesticide Mixture 236

DRE-GA09000236AL 1ml
Pesticide Mixture 6 600 µg/mL in Acetonitrile(‡)

acephate	chlorpyrifos
diazinon	dimethoate
ethoprophos (prophos)	malathion
methyl parathion	dibrom
phosmet	dichlorvos

Pesticide-Mix 168

DRE-L18000168TO 10ml
Pesticide-Mix 168 10 µg/mL in Toluene

Dicrotophos	Dimefox
Fenitrothion	Formothion
Methamidophos	Mevinphos
Parathion-ethyl	Parathion-methyl
Triazophos	Disulfoton
Dimethoate	Malaoxon
Iodofenphos	Paraoxon-ethyl
Omethoate	Pyrazophos
Phosalone	

Pesticide-Mix 195

DRE-LA18000195IO 1ml
Pesticide-Mix 195 10 µg/mL in Isooctane(‡)

Cypermethrin (technical) [10 µg/mL]	Deltamethrin [10 µg/mL]
Fenpropathrin [5 µg/mL]	Fenvalerate [10 µg/mL]
Permethrin [10 µg/mL]	

Pesticide Mixture 204

DRE-GH09000204AL 10x1ml
Pesticide Mixture 204 100 µg/mL in Acetonitrile(‡)(*)

acetamiprid	azoxystrobin
boscalid	buprofezin
chlorantraniliprole	chlorfluazuron
dicrotophos	Difenoconazole
dinotefuran	fenbuconazole
flonicamid	flubendiamide
hexaflumuron	imazail
bifenazate	bifenthrin
carbaryl	carbofuran
clothianidin	cyprodinil
Diflubenzuron	dimethoate
fenbutatin oxide	fipronil
fludioxonil	flufenoxuron

Pesticide-Mix 235

DRE-XA18000235AC 1ml
Pesticide-Mix 235 200 µg/mL in Acetone(‡)

Bromophos-ethyl	Bromophos-methyl
Chlorfenvinphos	Chlorpyrifos
Diazinon	Dichlorvos
Dimethoate	Disulfoton
Ethion	Fenthion
Malathion	Mevinphos
Paraoxon-ethyl	Parathion-ethyl
Parathion-methyl	

Pesticide Mixture 250

DRE-GS09000250AC 5x1ml
Pesticide Mixture 250 100 µg/mL in Acetone(‡)(*)

resmethrin	fenamiphos sulfoxide
fenamiphos sulfone	bifenthrin
tetramethrin	phosmet
etoxazole	fenamidone
triticonazole	lambda cyhalothrin
fenbuconazole	fenthion
tau-fluvalinate	

Pesticide Mixture 251

DRE-GS09000251AC 5x1ml
Pesticide Mixture 251 100 µg/mL in Acetone(±)(*)

flumioxazin	triadimefon
cyprodinil	fipronil
fenamiphos	flutolanil
tribufos	Systhane (TM)
carfentrazone-ethyl	hexazinon
fosthiazate	pendimethalin
allethrin	methidathion
fludioxonil	oxadiazon
buprofezin	ethion
piperyonyl butoxide	Propiconazol (mixture of isomers)

Pesticide Mixture 254

DRE-GS09000254AC 5x1ml
Pesticide Mixture 254 200 µg/mL in Acetonitrile(±)(*)

Lontrel (clopyralid)	chloramben
tricycopr	picloram
bentazon	2,4-dichlorophenoxyacetic acid
Dicamba	MCPP (Mecoprop)
MCPA acid	imazapyr
imazamox	imazaquin
imazethapyr	

Pesticide Mixture 301

DRE-GA09000301AL 10ml
Pesticide Mixture 301 100 µg/mL in Acetonitrile(±)

thiophanate methyl	3-hydroxycarbofuran
pyraclostrobin	abamectin
acetamiprid	amitraz
azoxystrobin	chlorantraniliprole
Imazalil	imidacloprid
methoxyfenozide	Systhane TM
piperyonyl butoxide	fenamidone
thiabendazole	

Pesticide Mixture 395

DRE-GS09000395AL 5x1ml
Pesticide Mixture 395 1-10 µg/mL in Acetonitrile(±)(*)

acephate [4 µg/mL]	aldicarb [4 µg/mL]
chlorpyrifos [2 µg/mL]	baythroid (four isomers) [10 µg/mL]
dimethoate [2 µg/mL]	ethoprophos (prophos) [2 µg/mL]
metalaxyl [2 µg/mL]	methiocarb [2 µg/mL]
oxamyl [10 µg/mL]	prallethrin [2 µg/mL]
azoxystrobin [2 µg/mL]	carbaryl [2 µg/mL]
dichlorvos [1 µg/mL]	diazinon [2 µg/mL]
kresoxim methyl [4 µg/mL]	malathion [2 µg/mL]
methomyl [4 µg/mL]	methyl parathion [2 µg/mL]
pyridaben [2 µg/mL]	cypermethrin (mix of isomers) [10 µg/mL]

Pesticide-Mix 1471

DRE-LA18001471CY 1ml
Pesticide-Mix 1471 10 µg/mL in Cyclohexane(±)

2,4'-DDT	4,4'-DDT
beta-HCH	cis-Chlordane (alpha)
Endrin	gamma-HCH
Heptachlor-exo-epoxide (cis-isomer B)	Hexachlorobenzene
Aldrin	alpha-HCH
Dicofol	Dieldrin
Heptachlor	Heptachlor-endo-epoxide (trans-isom. A)
oxy-Chlordane	trans-Chlordane (gamma)

Pesticide Mixture 253

DRE-GS09000253AC 5x1ml
Pesticide Mixture 253 100 µg/mL in Acetone(±)(*)

atrazine desisopropyl	desethyl atrazine
dimethoate	simazine
propyzamide (pronamide)	chlorothalonil
alachlor	malathion
endosulfan I	chlorfenapyr
iprodione	danilol
permethrin (mixture of isomers)	baythroid (mixture of isomers)
deltamethrin	phenothrin
Benfluralin (Benefin)	parathion
atrazine	g-BHC
metribuzin	methyl parathion
phorate	MGK-264 Mix of Isomers
fluzifop-P-butyl	endosulfan II
Guthion®	coumaphos
cypermethrin (mix of isomers)	fenvalerate (mixture of diastereoisomers)
prodiamine	

Pesticide Mixture 300

DRE-GA09000300AL 10ml
Pesticide Mixture 300 100 µg/mL in Acetonitrile(±)(*)

alpha-cypermethrin	allethrin
bifenthrin	baythroid (mixture of isomers)
danilol	deltamethrin
esfenvalerate	ethofenprox
fenvalerate (mixture of diastereoisomers)	lambda cyhalothrin
permethrin (mixture of isomers)	phenothrin

Pesticide-Mix 323

DRE-XA18000323IO 1ml
Pesticide-Mix 323 20-40 µg/mL in Isooctane(±)

2,4'-DDD [20 µg/mL]	2,4'-DDE [20 µg/mL]
4,4'-DDD [20 µg/mL]	4,4'-DDE [20 µg/mL]
alpha-Endosulfan [20 µg/mL]	alpha-HCH [20 µg/mL]
Chlorothalonil [40 µg/mL]	cis-Chlordane (alpha) [20 µg/mL]
Dieldrin [20 µg/mL]	Endosulfan-sulfate [20 µg/mL]
Heptachlor [20 µg/mL]	Heptachlor-exo-epoxide [20 µg/mL]
trans-Chlordane (gamma) [20 µg/mL]	2,4'-Dicofol [40 µg/mL]
2,4'-DDT [20 µg/mL]	Aldrin [20 µg/mL]
4,4'-DDT [20 µg/mL]	beta-HCH [20 µg/mL]
beta-Endosulfan [20 µg/mL]	Dicofol [20 µg/mL]
delta-HCH [20 µg/mL]	gamma-HCH [20 µg/mL]
Endrin [20 µg/mL]	oxy-Chlordane [20 µg/mL]
Hexachlorobenzene [20 µg/mL]	

Pesticide-Mix 1037

DRE-LA18001037CY 1ml
Pesticide-Mix 1037 10 µg/mL in Cyclohexane(±)

2,4'-DDD	2,4'-DDE
4,4'-DDE	4,4'-DDT
beta-HCH	cis-Chlordane (alpha)
Dieldrin	Endrin
Heptachlor-endo-epoxide (trans-isom. A)	Heptachlor-exo-epoxide (cis-isomer B)
trans-Chlordane (gamma)	trans-Nonachlor
2,4'-DDT	4,4'-DDD
Aldrin	alpha-HCH
cis-Nonachlor	delta-HCH
gamma-HCH	Heptachlor
Hexachlorobenzene	oxy-Chlordane

Pesticide-Mix 1584

DRE-LA18001584AL 1ml
Pesticide-Mix 1584 10 µg/mL in Acetonitrile

Atrazine	Atrazine-desethyl
Cyanazine	Dimethoate
Isoproturon	Linuron
Metribuzin	Pirimicarb
Propyzamide	Simazine
Triadimenol	Chloridazon
Atrazine-desisopropyl	Hexazinone
Diuron	Methabenzthiazuron
Metamitron	Propiconazole
Prochloraz	Terbutylazine-desethyl
Terbutylazine	

Pesticide-Mix 1598

DRE-LA18001598CY 1ml
Pesticide-Mix 1598 10 µg/mL in Cyclohexane

2,4'-DDT	4,4'-DDT
alpha-HCH	beta-Endosulfan
Dicofol	Dieldrin
Heptachlor	Heptachlor-endo-epoxide (trans-isom. A)
oxy-Chlordane	trans-Chlordane (gamma)
Aldrin	alpha-Endosulfan
beta-HCH	cis-Chlordane (alpha)
Endrin	gamma-HCH
Heptachlor-exo-epoxide (cis-isomer B)	Hexachlorobenzene

Pesticide Mixture 510

DRE-A50000510MB 1ml
Pesticide Mixture 510 1000 µg/mL in Methyl-tert-butyl ether(‡)

Parathion-methyl	Parathion-ethyl
Malathion	Dimethoate
Dichlorvos	Demeton (O+S)

Pesticides Mixture 596

DRE-A50000596AC 1ml
Pesticides Mixture 596 100 µg/mL in Acetone(‡)

carbofuran	3-hydroxycarbofuran
aldicarb	carbosulfan
aldicarb sulfone	aldicarb sulfoxide
methomyl	

Pesticide Mixture 679

DRE-A50000679AC 1ml
Pesticide Mixture 679 100 µg/mL in Acetone(‡)

pyrimethanil	mepanipyrim
Systhane TM	azoxystrobin

Pesticide Mixture B

DRE-GS09000052DI 4x1ml
Pesticide Mixture B in Dichloromethane(‡)

methomyl [1000 µg/mL]	Guthion® [1000 µg/mL]
chlorpyrifos [50 µg/mL]	dimethoate [50 µg/mL]
Ethalfuralin [50 µg/mL]	ethion [1000 µg/mL]
fenoxaprop-P-ethyl [250 µg/mL]	linuron [250 µg/mL]
terbufos [250 µg/mL]	napropamide [250 µg/mL]
bromacil [250 µg/mL]	carboxine [1000 µg/mL]
disulfoton [1000 µg/mL]	diuron [1000 µg/mL]
phorate [50 µg/mL]	thiamethoxam [250 µg/mL]
malathion [250 µg/mL]	pyridaben [250 µg/mL]
vinclozolin [50 µg/mL]	

Pesticide Mixture E

DRE-GS09000139TO 5x1ml
Pesticide Mixture E 100 µg/mL in Toluene(‡)(*)

acetochlor	alachlor
anilazine	atrazine
b-BHC	d-BHC
bifenthrin	Bioresmethrin
bromopropylate	butachlor
aldrin	ametryne
Benfluralin (Benefin)	a-BHC
g-BHC	bifenox
bitertanole (mixture of isomers)	boscalid
butralin	butylate

Pesticide-Mix 1612

DRE-LA18001612HE 1ml
Pesticide-Mix 1612 10 µg/mL in Hexane

4,4'-DDD	4,4'-DDE
2,4'-DDT	4,4'-DDT
alpha-HCH	beta-HCH
gamma-HCH	delta-HCH

Pesticides Mixture 12

DRE-B30000012AL 10ml
Pesticides Mixture 12 10 µg/mL in Acetonitrile

Ametryn	Atrazine
Atrazine-desethyl	Atrazine-desisopropyl
Cyanazine	Desmetryn
Hexazinone	Metribuzin
Prometryn	Propazine
Sebuthylazine	Simazine
Terbutylazine	Terbutryn

Pesticide Mixture 630

DRE-A50000630HE 1ml
Pesticide Mixture 630 100 µg/mL in Hexane(‡)

cypermethrin (mix of isomers)	lambda cyhalothrin
bifenthrin	fenpropathrin
baythroid (mixture of four isomers)	deltamethrin
fenvalerate (mixt. of diastereoisomers)	kelthane ® (dicofol)

Pesticide Mixture A

DRE-GS09000051DI 4x1ml
Pesticide Mixture A in Dichloromethane(‡)(*)

atrazine desisopropyl [1000 µg/mL]	chlorothalonil [200 µg/mL]
methoxychlor [1000 µg/mL]	metolachlor [200 µg/mL]
oxycarboxin [1000 µg/mL]	dieldrin [200 µg/mL]
g-BHC [200 µg/mL]	diazinon [200 µg/mL]
metalaxyl [200 µg/mL]	parathion [200 µg/mL]
bladex [1000 µg/mL]	desethyl atrazine [1000 µg/mL]
simazine [200 µg/mL]	aldicarb [4000 µg/mL]
aldrin [200 µg/mL]	a-BHC [200 µg/mL]
endosulfan I [200 µg/mL]	ethofumesate [200 µg/mL]
triallate [200 µg/mL]	

Pesticide Mixture C

DRE-GS09000053DI 4x1ml
Pesticide Mixture C in MeCl₂:Acetone:MeOH:MtBE 5:2:2:1(‡)(*)

brominal [200 µg/mL]	4-chloro-2-methylphenol [200 µg/mL]
2,4-dichlorophenol [200 µg/mL]	dichlorprop (2,4-DP) [200 µg/mL]
imazethapyr [1000 µg/mL]	MCPA acid [200 µg/mL]
picloram [200 µg/mL]	Triclopyr [200 µg/mL]
imazamox [200 µg/mL]	quinclorac [200 µg/mL]
2,4-DB [200 µg/mL]	Dicamba [200 µg/mL]
Lontrel (clopyralid) [1000 µg/mL]	imazamethabenz-methyl [1000 µg/mL]
diclofop-methyl [1000 µg/mL]	MCPP acid [200 µg/mL]
MCPB [1000 µg/mL]	fluroxypyr [200 µg/mL]
Hexaconazole [200 µg/mL]	2,4-D [200 µg/mL]
trifluralin [200 µg/mL]	

Pesticide Mixture D

DRE-GS09000054DI 4x1ml
Pesticide Mixture D in Dichloromethane(‡)

iprodione [1000 µg/mL]	clodinafop-propargyl [1000 µg/mL]
clodinafop (acid free) [1000 µg/mL]	fluzafop [200 µg/mL]
aminopyralid [1000 µg/mL]	Propiconazole (mixture of isomers) [1000 µg/mL]
quizalofop [1000 µg/mL]	bentazon [200 µg/mL]

Pesticide Mixture R

DRE-GS09000152AL 5x1ml
Pesticide Mixture R 100 µg/mL in Acetonitrile(±)(*)

pyridate	pyroxsulam
spinosad (Spinosyn A & D)	sulfentrazone
terbutol	thiabenzazole
thifensulfuron-methyl	thiobencarb
thiofanox-sulfone	thiophanate methyl
rimsulfuron	sethoxydim
tebufenozide	tembotrione
thiacloprid	thiamethoxam
thiodicarb	thiofanox
tolyfluanid	triasulfuron

Pesticide Surrogate Mixture 541

DRE-A50000541TH 1ml
Pesticide Surrogate Mixture 541 2000 µg/mL in Toluene:Hexane(±)

dibutyl chlorendate	decachlorobiphenyl (BZ# 209)
2,4,5,6-tetrachloro-m-xylene	

Pesticides Mixture 497 for HJ 961-2018, HJ 1026-2019

DRE-A50000497ME 1ml
HJ 961-2018, HJ 1026-2019 Pesticides Mixture 497 100 µg/mL in Methanol(±)

Oxamyl	Methomyl
Dioxacarb	Aldicarb
Bendiocarb	Carbofuran
Propoxur	Carbaryl
Ethiofencarb	Pirimicarb
Isoprocarb	Fenobucarb
Methiocarb	Promecarb
Alanycarb	

Rodenticides Mixture 248

DRE-GS09000248AL 5x1ml
Rodenticides Mixture 248 100 µg/mL in Acetonitrile(±)

bromadiolon	diphacinone
chlorophacinone	Brodifacoum
difenacoum	warfarin

Smart Solutions™ v400 GC PestiMix Kit 1

DRE-K50000291IT 1kit (5 ea - 1x1ml)
Smart Solutions™ v400 GC PestiMix Kit 1(±)(*)

DRE-A50000292I	GC PestiMix 1 10 µg/mL in Isooctane:Toluene (50:50)
DRE-A50000293I	GC PestiMix 2 10 µg/mL in Isooctane:Toluene (50:50)
DRE-A50000294I	GC PestiMix 3 10 µg/mL in Isooctane:Toluene (50:50)
DRE-A50000295I	GC PestiMix 4 10 µg/mL in Isooctane:Toluene (50:50)
DRE-A50000296I	GC PestiMix 5 10 µg/mL in Isooctane:Toluene (50:50)

Pesticide Mixture S

DRE-GS09000153AL 5x1ml
Pesticide Mixture S 100 µg/mL in Acetonitrile(±)(*)

tribenuron-methyl	trichlorfon
tridemorph (mix of isomers)	trifloxystrobin
triflumizole	XMC (3,5-xylyl methylcarbamate)

Pesticide Surrogate Mixture 489 for HJ 699-2014, HJ 904-2017, HJ 901-2017, HJ 903-2017, HJ 743-2015

DRE-A50000489TH 1ml
HJ 699-2014, HJ 904-2017, HJ 901-2017, HJ 903-2017,
HJ 743-2015 Pesticide Surrogate Mixture 489
1000 µg/mL in Toluene:Hexane(±)

PCB 209 (Decachlorobiphenyl)	2,4,5,6-Tetrachloro-m-xylene
------------------------------	------------------------------

Pesticides Mixture 484 for HJ 1052-2019

DRE-A50000484AL 1ml
HJ 1052-2019 Pesticides Mixture 484 100 µg/mL in Acetonitrile(±)

Simazine	Atraton
Simetryn	Atrazine
Secbumeton	Prometon
Ametryn	Propazine
Terbutylazine	Prometryn
Terbutryn	

Pesticides Surrogate Standard Spiking Solution

DRE-GA09000914AC 1ml
Pesticides Surrogate Standard Spiking Solution 200 µg/mL in Acetone(±)

decachlorobiphenyl (BZ# 209)	2,4,5,6-tetrachloro-m-xylene
------------------------------	------------------------------

Pyridine Pesticide Mixture 675

DRE-A50000675ME 1ml
Pyridine Pesticide Mixture 675 500 µg/mL in Methanol(±)

fluroxypyr	dithiopyr
thiazopyr	

Single Column Analytes Mix 4

DRE-XA08080400TH 1ml
Single Column Analytes Mix 4 100 µg/mL in Toluene/Hexane
DRE-YA08080400TH 1ml
Single Column Analytes Mix 4 2000 µg/mL in Toluene/Hexane

4,4'-DDD	4,4'-DDE
alpha-Endosulfan	alpha-HCH
cis-Chlordane (alpha)	delta-HCH
Endrin	Endrin-aldehyde
Heptachlor	Heptachlor-endo-epoxide (trans-isom. A)
4,4'-DDT	Aldrin
beta-Endosulfan	beta-HCH
Dieldrin	Endosulfan-sulfate
Endrin-ketone	gamma-HCH
trans-Chlordane (gamma)	Methoxychlor

GC PestiMix 1

DRE-A50000292IT 1ml
GC PestiMix 1 10 µg/mL in Isooctane:Toluene (50:50)(±)

Acephate	Atraton
Benalaxyl	Bifenthrin
Butralin	Cadusafofos
Chlorfenvinphos	Chlorobenzilate
Chlorpyrifos	Chlorpyrifos-methyl
Cyprodinil	4,4'-DDD
Atrazine	Atrazine-desisopropyl
Bitertanol	Bromophos-ethyl
Cis-Chlordane (Alpha Isomer)	Chlorfenapyr
Chloroneb	Chlorpropham
Chlorthal-dimethyl	Cyproconazole
Deltamethrin	Demeton-S-methyl

(continued on next page)

GC PestiMix 2

(continued from previous page)

DRE-A50000293IT 1ml
GC PestiMix 2 10 µg/mL in Isooctane:Toluene (50:50)(±)

Acetochlor	Aclonifen
Azinphos-ethyl	Azinphos-methyl
Biphenyl	Bromacil
Butylate	Butylhydroxytoluene
Carboxin	Chlormephos
Cypermethrin	2,4-D 2-ethylhexyl ester
Dichlorvos	Dicofol
Dimethomorph	2,4-Dinitrotoluene
Disulfoton-sulfone	Edifenphos
Famphur	Fenchlorphos
Fensulfothion	Fenthion-sulfone
beta-HCH	Heptachlor-endo-epoxide (isom. A)
Indoxacarb	Isodrin
Kresoxim-methyl	Malathion
Metolachlor	Monocrotophos
trans-Nonachlor	Norflurazon
Pebulate	Penconazole
Phosphamidon	Pirimiphos-methyl
Promecarb	Propyzamide
Tebupirimfos	Terbacil
Thionazin	Triadimefon
Tricyclazole	O,O,O-Triethylphosphorothioate
Acrinathrin	Alachlor
Bendiocarb	Bifenox
1-Bromo-2-nitrobenzene	Bromophos-methyl
Carbaryl	Carbofuran
Chlorothalonil	Coumaphos
4,4'-DDT	Demeton (O+S)
Diethyltoluamide (DEET)	Dimethenamid
2,6-Dinitrotoluene	9,10-Diphenylanthracene
Endrin	EPTC
Fenitrothion	(±)-Fenprothrin
Fluridone	Fluroxypr-1-methylheptylester
Hexachlorocyclopentadiene	Imazalil
Isofenphos-methyl	Isoprothiolane
Methamidophos	Methidathion
Naled	Nitrofen
Omethoate	Parathion-ethyl
Phosalone	Phosmet
Prochloraz	Profuralin
Sulprofos	tau-Fluvalinate
Terbumeton	Tetrachlorvinphos
Triallate	Tributylphosphate
Triphenylphosphate	Vernolate

Smart Solutions™ v700 LC PestiMix Kit 2

DRE-K50000087 1set (5ea x 1x1ml)
Smart Solutions™ v700 LC PestiMix Kit 2(±)(*)

DRE-A50000082CMLC PestiMix 1
5 µg/mL in Acetone:Acetonitrile:Methanol (50:49:1)

DRE-A50000083AALC PestiMix 2
5 µg/mL in Acetonitrile:Acetone (94.5:5.5)

DRE-A50000084AALC PestiMix 3
5 µg/mL in Acetonitrile:Acetone (69:21)

DRE-A50000085AALC PestiMix 4
5 µg/mL in Acetonitrile:Acetone (72:13.5)

DRE-A50000086ALLC PestiMix 5
5 µg/mL in Acetonitrile

GC PestiMix 4

DRE-A50000295IT 1ml
GC PestiMix 4 10 µg/mL in Isooctane:Toluene (50:50)(±)

Acequinocyl	Adipic Acid Bis(2-ethylhexyl) Ester
Aminocarb	Bentazone
Butachlor	Butamifos
Chlozolinate	Cyanophos
Dibutyl phthalate	3,4-Dichloroaniline
Dieldrin	Diffufenican
2,6-Dimethylphenol	3,4-Dimethylphenol
Endrin-ketone	Epoxiconazole
Ethofumesate	Fenamiphos
Fenthion-sulfoxide	Fipronil

(continued on next page)

Diazinon	1,2-Dibromo-3-chloropropane
Dichlofluanid	Difenoconazole
(EZ)-Diniconazole	Diphenylamine
Etofenprox	Etozazole
Fenbuconazole	Fenoxaprop-P-ethyl
Fluquinconazole	Flusilazole
Heptenophos	Hexachlorobenzene
Iprobenfos	Isoprocarb
Mevinphos	Mirex
Nuarimol	Oxadiazon
2-Phenylphenol	Phorate
Prometryn	Propachlor
Prosulfocarb	Pyrazophos
Pyrimethanil	Pyriproxyfen
Tebufenpyrad	Tefluthrin
2,4,5,6-Tetrachloro-m-xylene	Tetraconazole
Triadimenol	Triazophos
Dichlobenil	Dichlofenthion
Dimethoate	Diniconazole
Disulfoton-sulfoxide	Ethoprophos
Etridiazole	Fenazaquin
Fenvalerate	Flucythrinate
alpha-HCH	Heptachlor
Hexaconazole	Hexazinone
Mefenpyr-diethyl	Methacrifos
Molinate	Napropamide
Paraoxon-methyl	Phenthoate
Phosfolan	Profenofos
Propazine	Propiconazole
Pyridaphenthion	Pyrifenoxy
Simazine	Sulfotep
Terbufos	2,3,5,6-Tetrachloronitrobenzene
Tetradifon	Thiometon
Trifluralin	Vinclozolin

GC PestiMix 3

DRE-A50000294IT 1ml
GC PestiMix 3 10 µg/mL in Isooctane:Toluene (50:50)(±)

Acetamiprid	Ametoctradin
Benfluralin	Boscalid
Captan	Carbetamide
Chlorthiophos	Clodinafop-propargyl ester
Cycloate	Cyfluthrin
2,4'-DDD	2,4'-DDE
Diclobutrazol	Diclofop methyl
Diphenamid	Disulfoton
Endrin-aldehyde	Esfenvalerate
Famoxadone	Fenoxycarb
Fluchloralin	Flufenacet
Folpet	Fonofos
Isofenphos	Mepronil
Methiocarb	Methoxychlor
cis-Nonachlor	Oxamyl
Parathion-methyl	PCB 209
Permethrin	Perthane
Pirimicarb	Pirimiphos-ethyl
Propamocarb free base	Propanil
Pyraclifos	Pyraclostrobin
Simetryn	S-Metolachlor
Thiazopyr	Trichloronate
Atrazine-desethyl	Azoxystrobin
Bromuconazole	Captafol
Carbophenothion	Carfentrazone-ethyl
Clomazone	Cyanazine
lambda-Cyhalothrin	2,4-D-iso-propyl ester
4,4'-DDE	Diallate
Dicloran	Dioxacarb
Endosulfan-sulfate	alpha-Endosulfan
Ethalfuralin	Etrifmos
Fenpropidin	Fluazifop-P-butyl
Flumorph	Flutolanil
gamma-HCH	Ipconazole
Metalaxyl-M	Metazachlor
Metolcarb	Myclobutanil
oxy-Chlordane	Paraoxon-ethyl
Pendimethalin	Pentachlorobenzene
Phorate-sulfone	Piperonyl butoxide
Pretilachlor	Prometon
Propham	Propoxur
Pyridaben	Quinalphos
Spiromesifen	Terbutylazine
Trifloxystrobin	

(continued from previous page)

delta-HCH	Heptachlor-exo-epoxide (isom. B)
Isazofos	Isobenzan (Telodrin)
Methomyl	Metribuzin
Oxyfluorfen	Pentachloroanisole
Procymidone	Prothiophos
Tebutam	Terbutryn
1,2,4,5-Tetrachlorobenzene	2,3,4,6-Tetrachlorophenol
Tolyfluandil	Tralomethrin
Aldrin	Ametryn
Bromopropylate	Bromoxynil
trans-Chlordane (Gamma Isomer)	Chlordene
Cymiazole	2,4'-DDT
2,4-Dichlorophenol	Dicrotophos
Dimethachlor	Dimethipin
Dinoterb	beta-Endosulfan
Ethiofencarb	Ethion
Fenpropathrin	Fenpropimorph
Flurtamone	Fosthiazate
Hexachlorobutadiene	Hexachloroethane
Isocarbofos	Metaxyl
Nitrothal-isopropyl	Norflurazon-desmethyl
Picoxystrobin	Prallethrin
Quintozene	Tebuconazole
1,2,3,4-Tetrachlorobenzene	1,2,3,5-Tetrachlorobenzene
Thiobencarb	Tolclofos-methyl
2,4,6-Trichlorophenol	Triticonazole

PestiMix 1

DRE-A50000082CM 1ml

LC PestiMix 1 5 µg/mL in Acetone:Acetonitrile:Methanol (50:49:1)(‡)(*)

Acetamidiprid	Acibenzolar-S-methyl
Aldicarb-sulfoxide	Allidochlor
Atraton	Atrazine
Benalaxyl	Bendiocarb
Benzoximate	Bitertanol
Butafenacil	Butocarboxim
Carbaryl	Carbendazim
Carbofuran-3-hydroxy	Carboxin
Chlorbufam	Chlorfluazuron
Chloroxuron	Chlorsulfuron
Clodinafop-propargyl ester	Cyanazine
Aldicarb	Aldicarb-sulfone
Ametoctradin	Ametryn
Atrazine-desethyl	Azoxystrobin
Bensulfuron-methyl	Benthiavicalcarb-isopropyl
Boscalid	Buprofezin
Butoxycarboxim	Buturon
Carbetamide	Carbofuran
Chlorantraniliprole	Chlorbromuron
Chloridazon	Chlormequat chloride
Cinosulfuron	Clethodim
Cyazofamid	Cycluron
Cymoxanil	Cyprodinil
Dichlormid	Diethofencarb
Dimethomorph	Dinotefuran
Diuron	DMST
Ethiprole	Ethirimol
Etoxazole	Famoxadone
Fenobucarb	Fenoxaprop-P-ethyl (R-enantiomer)
Fenthion-sulfone	Fenuron
Fluazifop-P-butyl	Flubendiamide
Fluoxastrobin	Flusilazole
Formetanate hydrochloride	Fuberidazole
Hexythiazox	Imazalil
Imazosulfuron	Imidacloprid
Iprodione	Iprovalicarb
Isoxaflutole	Lenacil
Malathion	Mandipropamid
Mepanipyrim	Mephosfolan
Mesosulfuron-methyl	Methabenzthiazuron
Methidathion	Methiocarb
Metobromuron	Metolcarb
Mevinphos	Mexacarbate
Monuron	Moxidectin
Nuarimol	Ofurace
Paraoxon-ethyl	Paraoxon-methyl
Phorate	Phosalone
Phoxim	Piperonyl butoxide

(continued on next page)

GC PestiMix 5

DRE-A50000296IT 1ml

GC PestiMix 5 10 µg/mL in Isooctane:Toluene (50:50)(‡)(*)

Acibenzolar-S-methyl	Allethrin
Azaconazole	Benzoylprop-ethyl
Bupirimate	Chlorbencide
Chlorfensan	1-Chloro-3-nitrobenzene
gamma-Cyhalothrin	beta-Cypermethrin
4,4'-DDMU	Desmedipham
4,4'-Dichlorobenzophenone	1,3-Dimethyl-2-nitrobenzene
epsilon-HCH	Ethoxyquin
Fenarimol	Fenchlorphos-oxon
Fenson	Fenthion
Fludioxonil	Fluometuron
Halfenprox	(S)-Indoxacarb
Irgarol 1051	Isopropalin
Linuron	Malaoxon
2,4'-Methoxychlor	4,4'-Methoxychlor olefin
Nitralin	Oxycarboxin
Pentachlorobenzonitrile	cis-Permethrin
Phthalic acid bis-2-ethylhexyl ester	Picolinafen
Propetamphos	Propisochlor
Quinoxifen	Resmethrin
Spiroxamine	Tebuthiuron
cis-1,2,3,6-Tetrahydrophthalimide	Tetramethrin
Transfluthrin	Triflumizole
Allidochlor	Anthraquinone
Bromfenvinfos-methyl	Bromfenvinphos
Chlordecone	Chlordimeform free base
Chloropropylate	Cloquintocet-1-methylhexyl ester
Dazomet	DCAA Methyl Ester
Dibutyl chlorendate	Dicaphon
Endosulfan-ether	EPN
Fenamidone	Fenamiphos-sulfone
Fenhexamid	Fenoxaprop-ethyl
Flonicamid	Fluazifop-butyl
Flutriafol	Form-2',4'-Xylidide
Iodofenphos	Iprodione
Lenacil	Leptophos
Mecarbam	Methiocarb sulfone
Methyl-pentachlorophenyl sulfide	MGK 264
Pacloutrazol	Pentachloroaniline
trans-Permethrin	Phenothrin
Prodiamine	Propargite
Prothoate	(Z)-Pyrinobac-methyl
S 421	Spirodiclofen
O,O-TEPP	2,3,5,6-Tetrachloroaniline
Thiamethoxam	Tiocarbazil
Zoxamide	

PestiMix 2

DRE-A50000083AA 1ml

LC PestiMix 2 5 µg/mL in Acetonitrile:Acetone (94.5:5.5)(‡)(*)

Abamectin	Acephate
Allethrin	Azaconazole
Azinphos-methyl	Beflubutamid
Benzoylprop-ethyl	Bifenthrin
Bupirimate	Butachlor
Carfentrazone-ethyl	Chinomethionate
Chlorthiophos	Chromafenozide
Clofentazine	Clomazone
Cyanofenphos	Cyanophos
Cypermethrin (technical)	Cyphenothrin
Demeton-S-methyl	Demeton-S-methyl sulfone
Diazinon	Dichlofenthion
Dicrotophos	Difenoconazole
Dimethenamid	Dimethoate
Diniconazole	Dioxathion
Disulfoton-sulfoxide	Ditalimfos
Empenthrin	EPN
Etaconazole	Ethiofencarb-sulfone
2-Ethylhexyl diphenyl phosphate	Etrifmos
Fenamiphos-sulfone	Fenamiphos-sulfoxide
Fenbuconazole	Fenchlorazol-ethyl
Fenoxaprop	Fenpropimorph
Fensulfthion-oxon-sulfone	Fensulfthion-sulfone
Fipronil	Flamprop-isopropyl
Fluazinam	Flucythrinate

(continued on next page)

Pirimicarb-desmethyl	Pirimicarb-desmethyl-formamido	Fluquinconazole	Flurenol-butyl ester
Profenofos	Promecarb	Fluroxypyr-1-methylheptylester	Flurprimidol
Propaquizafop	Propargite	Fosthiazate	Furathiocarb
Propham	Propulfocarb	Haloxyfop-ethoxyethyl	Haloxyfop-methyl
Pyracarbolid	Pyrazophos	Aclonifen	Alachlor
Quinalphos	Quinchlorac	Azamethiphos	Azinphos-ethyl
Sebuthylazine	Siduron	Benoxacor	Bentazone
Spinetoram	Spinosad	Bixafen	Bromfenvinphos
Spiroxamine	Sulfotep	Cadusafos	Carbophenothion
Tebufenpyrad	Tebuthiuron	Chlorpyrifos	Chlorpyrifos-methyl
Terbufos	Terbumeton	Cinidon-ethyl	Climbazole
Terbutryn	Tetrachlorvinphos	Cloquintocet-1-methylhexyl ester	Coumaphos
Thiamethoxam	Thidiazuron	Cycloxydim	lambda-Cyhalothrin
Thiodicarb	Thiofanox	Cyproconazole	Cyprofuram
Triadimenol	Triasulfuron	Demeton-S-methyl sulfoxide	Dialifos
Tricyclazole	Tridemorph	Dichlorvos	Diclobutrazol
Triflumuron	Vamidothion	Diflubenzuron	Dimethachlor
Cyromazine	Desmedipham	Dimethylvinphos	Dimoxystrobin
Diethyltoluamide (DEET)	Diflufenican	Diphenamid	Disulfoton-sulfone
Dioxacarb	Dipropetryn	Doramectin	Edifenphos
Ethiofencarb	Ethiofencarb-sulfoxide	Epoxiconazole	Eprinomectin
Ethoxysulfuron	Etofenprox	Ethion	Ethoprophos
Fenamidone	Fenitrothion	Famphur	Fenamiphos
Fenoxycarb	(E)-Fenpyroximate	Fenarimol	Fenazaquin
Flonicamid	Florasulam	Fenfuram	Fenhexamid
Flufenoxuron	Flumioxazin	Fensulfothion	Fensulfothion-oxon
Flutolanil	Forchlorfenuron	Fenthion	Fenthion-sulfoxide
Furalaxyl	Halofenozide	Flamprop-methyl	Fluacrypyrim
Imazamethabenz-methyl	Imazamox	Flufenacet	Fluotrimazole
(S)-Indoxacarb	Iodosulfuron-methyl sodium	Flurenol-methyl ester	Flurochloridone
Isoprocarb	Isoxaben	Flutriafol	tau-Fluvalinate
Linuron	Malaoxon	Furmecyclox	Haloxyfop (free acid)
Mecarbam	Mefenacet	Heptenophos	Hexaconazole
Mepiquat chloride	Mepronil	Hexaflumuron	Hexazinone
Methacrifos	Methamidophos	Isazofos	Isocarbofos
Methoprotryne	Methoxyfenozide	Isxadifen-ethyl	Ivermectin
Metoxuron	Metsulfuron-methyl	Mefenpyr-diethyl	Metaflumizone
Monocrotophos	Monolinuron	Metazachlor	Metconazole
Myclobutanil	Novaluron	Methomyl	S-Metolachlor
Omethoate	Orbencarb	Molinate	Napropamide
Pencycuron	Phenthoate	Oryzalin	Oxadiazyl
Phosmet	Phosphamidon	Oxamyl	Oxamyl-oxime
Piperophos	Pirimicarb	Pendimethalin	Phenmedipham
Pirimiphos-ethyl	Pirimiphos-methyl	Picolinafen	Picoxystrobin
Prometon	Prometryn	Prochloraz	Prochloraz
Propazine	Propetamphos	Propiconazole	Propoxur
Prosulfuron	Prothiophos	Prothioconazole-desthio	Pyraclostrobin
Pyridaphenthion	Pyriproxyfen	Pyridalyl	Pyrimethanil
Quizalofop-ethyl	Rotenone	(5EZ)-Selamectin	Spirotetramat
Simazine	Simetryn	Sulfosulfuron	Tebutam
Spirodiclofen	Spiromesifen	Terbufos-sulfoxide	Tetraconazole
Tebuconazole	Tebufenozide	Thiofanox-sulfoxide	Tiocarbazil
Temephos	Tepraloxymid	Triadimefon	Triallate
Terbuthylazine	Terbuthylazine-desethyl	Triforine	2,3,5-Trimethacarb
Thiabendazole	Thiacloprid	Imibenconazole	Ipconazole
Thifensulfuron-methyl	Thiobencarb	Isofenphos	Isoprothiolane
Thionazin	Tolclofos-methyl	Kresoxim-methyl	Lufenuron
Triazophos	Trichlorphon	Metalaxyl	Metamitron
Trietazine	Trifloxystrobin	Methiocarb sulfone	Methiocarb sulfoxide
XMC (3,5-xylil methylcarbamate)	Zoxamide	Metrafenone	Metribuzin

PestiMix 3

DRE-A5000084AA 1ml

LC PestiMix 3 5 µg/mL in Acetonitrile:Acetone (69:21)(±)(*)

Aconitine	Amicarbazone
Asulam	Atrazine-desisopropyl
Benfluralin	Benfuresate
Bensulide	Benzthiazuron
Bromophos-methyl	Bromuconazole
Butoxydim	Butylate
Chlorfenvinphos	Chlorimuron-ethyl
Chlorthiamid	Cinmethylin
Cloransulam-methyl	Clothianidin
Crotoxyphos	Crufomate
Cycloate	Cycloheximide
Cymiazole	Cyprazine
Daimuron (Dymron)	Daminozide
Desmetryn	Diazoxon

(continued on next page)

(continued from previous page)

2,4-Dimethylaniline	Dimetilan
Diphenyl Phosphate	Dithiopyr
Emamectin	EPTC
Ethametsulfuron-methyl	Ethephon
Ethylene thiourea	Etobenzanid
Fenoxanil	Fenpiclonil
Fenthoxon Sulfone	Fenthoxon Sulfoxide
Flazasulfuron	Fluazifop (free acid)
Flucarbazone-sodium	Fluchloralin
Flumetralin	Flumorph
Fluridone	Flurtamone
Haloxypop-R-methyl	Hymexazol
Imidacloprid-olefin	Inabenfide
Isofenphos-methyl	Isofenphos-oxon
Isouron	Isoxathion
Levofloxacin	Mefluidide
Metosulam	Morpholine
Nitralin	Nitroguanidine
Octhilinone	Oxasulfuron
Parathion-methyl	Pebulate
Pentachlor	Penthiopyrad
Phorate Oxon Sulfone	Phorate Oxon Sulfoxide
Phosfolan	Pirimiphos-methyl-N-desethyl
Procyomidone	Profoxydim
Pyraclafos	Pyrazolynate
Pyrifenoxy	Pyroxulam
Quinoline	Quizalofop free acid
Roxithromycin	Saflufenacil monohydrate
Simeconazole	Simeton
Sulfamerazine	Sulfuramid
Tembotrione	O,O-TEPP
Thiazafurion	Thiazopyr
Tiadnil	Tolyfluanid
Tribenuron-methyl (technical)	Tri-o-cresyl phosphate
Triflusaluron-methyl	3,4,5-Trimethacarb
Vernolate	Aspon
Anilofos	Benazolin
Barban	Benomyl
Benodanil	BDMC
Bromacil	Butocarboxim sulfoxide
Bufencarb	Chlorfenapyr
Carpropamid	Chlorpyrifos-oxon
Chlorobenzuron	Clodinafop free acid
Clarithromycin	Crimidine
Coumaphos-oxon	Cyclanilide
Cumyluron	Cyflufenamid
Cyclosulfamuron	Cythioate
Cyrosulfamide	Demeton-O
Dazomet	Dimepiperate
Dimefox	Dinitramine
(EZ)-Diniconazole	Dodemorph
Dodecylguanidinium Acetate	Esprocarb
Erythromycin A	Ethoxyquin
Ethidimuron	Fenothiocarb
Fenchlorphos-oxon	Fenthoxon
Fenthion-ethyl	Fipronil Sulfide
Fentin-acetate	Fluazuron
Fluazifop-butyl	Fludioxonil
Flucycloxuron	Flupyradifurone
Fluopicolide	Furilazole
Fonofos	Imazethapyr
Imazaquin	Iprobenfos
Iodofenphos	Isopropalin
Isonoruron	Leptophos
Karbutilate	Metolachlor
Metalaxyl-M	Naled
Naftalofos	Norflurazon-desmethyl
Norflurazon	Parathion-ethyl
Oxaziclomefone	Penoxsulam
Penfluron	N-Phenylurea
Pethoxamid	Phorate-sulfoxide
Phorate-sulfone	Procyazine
Probenazole	Propisochlor
Propamocarb hydrochloride	Pyributicarb
Pyrazosulfuron-ethyl	Quinoclamine
Quinine sulfate	Rimsulfuron
Reserpine	Sethoxydim
Sebutylazine-desethyl	Sulfallate
Sulcotrione	Tebupirimfos
Sulfometuron-methyl	Terbucarb

(continued on next page)

GC PestiMix 5

DRE-A50000296IT	1ml
GC PestiMix 5	10 µg/mL in Isooctane:Toluene (50:50)(‡)(*)
Akton	Alloxydim-sodium
Aziprotryne	Bentazone-methyl
Bioresmethrin	Cafenstrole
Clomeprop	Clopyralid
alpha-Cypermethrin	Deoxynivalenol
1-(3,4-Dichlorophenyl)-3-methyl urea	Diclocymet
2,6-Diethylaniline	Dimethametryn
N,N'-Diphenylurea	Ethaboxam
Fenazox	Fenclorim
Fipronil-desulfinyl	Flamprop-M-isopropyl
Fluopyram	Fluorodifen
Fluroxypyr	Fluthiacet-methyl
Imazapic	Imiprothrin
Irgarol 1051	Isomethiozin
Methfuroxam	Methomyl-oxime
Z-Mevinphos (trans-butenoic acid)	MGK 264
(E)-Nitenpyram	Nitrapyrin
Oxycarboxin	Oxyfluorfen
Phospholan-methyl	Piperalin
Pyrazoxyfen	Pyrethrins
Schradan	Secbumeton
Sulfoxaflor	Sulfosulfoxide
Terbutylazine-2-hydroxy	Thiocyclam hydrogenoxalate
Triaziflam	Triazoxide
Trinexapac-ethyl	Tris(1,3-dichloroisopropyl) phosphate
Ancymidol	Azadirachtin A
1,2-Benzisothiazol-3(2H)-one	S-Bioallethrin
Carbofuran-3-keto	6-Chloro-4-hydroxy-3-phenyl-pyridazin
Cyantraniliprole	Cycloprothrin
3,4-Dichloroaniline	2,6-Dichlorobenzamide
Diclosulam	Diethatyl-ethyl
Dimethirimol	Diphenylamine
Ethiozin	Ethyachlozate
Fenpropathrin	Ferimzone
Flumetsulam	Flumiclorac-pentyl
Fluoroglycofen-ethyl	Flupyrasulfuron-methyl sodium
Fomesafen	S-Hydroprene
Indanofan	Lactofen
Ketoprofen	(Z)-Metaminostrobin
(E)-Metaminostrobin	Nereistoxin oxalate
Naptalam	Norastrobilin
Noruron	Phosmetoxon
Pentachlorobenzonitrile	Propazine-2-hydroxy
Prohydrojasmon	Quinmerac
Pyroquilon	Silthiofarn
Silthiofarn	Sulfamonomethoxine
Terbumeton-desethyl	Terbutaline Sulfate
Thiosultap-sodium	Tiamulin fumarate
Tridiphane	Trifloxystrobin (free acid)
Vamidothion-sulfoxide	

PestiMix 5

DRE-A50000086AL	1ml
LC PestiMix 5	5 µg/mL in Acetonitrile(‡)(*)
Amidosulfuron	Aminocarb
Chlordimeform	Chlormephos
Danitol	Demeton S
Dimefuron	Dithianon
Fenpropidin	Fluometuron
Isocarbamide	Isoproturon
Primisulfuron-methyl	Propamocarb
Pyridate	Quinoxifen
Thiencarbazone-methyl	Thiophanate methyl
Anilazine	Azimsulfuron
Chlorotoluron	Cyhalofop-butyl
Diafenthiuron	Difenoxurone
Ethofumesate	Fenazaflor
Foramsulfuron	Halosulfuron-methyl
Mesotrione	Nicosulfuron
Prothioconazole	Pymetrozine
Spirotetramat-enol	Spirotetramat-enol-glucoside
Tritosulfuron	

(continued from previous page)

Thifluzamide	Thenylchlor
Topramezone	Thiophanate-ethyl
Trifloxysulfuron Sodium Salt	Triazamate
Tris-(2-chloroethyl)phosphate	Trifluralin
	Uniconazole

Triazine Pesticides Mixture 926

DRE-GA09000926AC 1ml

Triazine Pesticides Mixture 926 100 µg/mL in Acetone(±)

ametryne	atraton
atrazine	prometon
prometryn	propazine
secbumeton	simetryn
simazine	terbuthylazine
terbutryne	

UCMR 4 Method 525.3

DRE-GS09000487ME 5x1ml

UCMR 4 Method 525.3, 10000 X MRL in Methanol(±)

a-BHC [100 µg/mL]	chlorpyrifos [300 µg/mL]
dimethipin [2000 µg/mL]	ethoprophos (prophos) [300 µg/mL]
oxyfluorfen [500 µg/mL]	tebuconazol (Folicur) [2000 µg/mL]
permethrin (mixture of isomers) [400 µg/mL]	tribufos [700 µg/mL]
profenofos [3000 µg/mL]	

Triazine & Urea Pesticide Mixture 447

DRE-A50000447AL 1ml

Triazine & Urea Pesticide Mixture 447 100 µg/mL in Acetonitrile(±)(*)

Atrazine	Atrazine-desethyl
Chloridazon	Metoxuron
Simazine	Cyanazine
Chlortoluron	Monolinuron
Metobromuron	Metazachlor
Terbuthylazine	Linuron
Chlorpropham	Terbutryn
Ethidimuron	Metamitron
Atrazine-desisopropyl	Bromacil
Carbetamide	Methabenzthiazuron
Terbuthylazine-desethyl	Isoproturon
Diuron	Dimefuron
Propazine	Prometryn
Chloroxuron	Ethofumesate
Metolachlor	

Washington Pesticide Mixture 1

DRE-A50000143AL 1ml

Washington Pesticide Mixture 1 1-20 µg/mL in Acetonitrile(±)(*)

DRE-V50000143AL 5x1ml

Washington Pesticide Mixture 1 1-20 µg/mL in Acetonitrile(±)(*)

Abamectin [5 µg/mL]	Acephate [4 µg/mL]
Aldicarb [4 µg/mL]	Azoxystrobin [2 µg/mL]
Boscalid [4 µg/mL]	Carbaryl [2 µg/mL]
Chlorfenapyr [10 µg/mL]	Chlormequat chloride [1 µg/mL]
Cyfluthrin [10 µg/mL]	Cypermethrin (technical) [10 µg/mL]
Dichlorvos [1 µg/mL]	Dimethoate [2 µg/mL]
Etoazole [2 µg/mL]	Fenoxycarb [2 µg/mL]
Flonicamid [10 µg/mL]	Fludioxonil [4 µg/mL]
Imidacloprid [4 µg/mL]	Kresoxim-methyl [4 µg/mL]
Methiocarb [2 µg/mL]	Methomyl [4 µg/mL]
Naled [5 µg/mL]	Oxamyl [10 µg/mL]
Permethrin [2 µg/mL]	Phosmet [2 µg/mL]
Propiconazole [4 µg/mL]	Propoxur [2 µg/mL]
Spinosad [2 µg/mL]	Spiromesifen [2 µg/mL]
Tebuconazole [4 µg/mL]	Thiacloprid [2 µg/mL]
Uniconazole [1 µg/mL]	Acetamiprid [2 µg/mL]
Acequinocyl [20 µg/mL]	Bifenthrin [2 µg/mL]
Bifenazate [2 µg/mL]	Chlorantraniliprole [2 µg/mL]
Carbofuran [2 µg/mL]	Clofentezine [2 µg/mL]
Chlorpyrifos [2 µg/mL]	Diazinon [2 µg/mL]
Daminozide [10 µg/mL]	Etofenprox [4 µg/mL]
Ethoprophos [2 µg/mL]	Fipronil [4 µg/mL]
(E)-Fenpyroximate [4 µg/mL]	Imazalil [2 µg/mL]
Hexythiazox [10 µg/mL]	Metalaxyl [2 µg/mL]
Malathion [2 µg/mL]	Myclobutanil [2 µg/mL]
MGK 264 [2 µg/mL]	Parathion-methyl [2 µg/mL]
Paclobutrazol [4 µg/mL]	Prallethrin [2 µg/mL]
Piperonyl butoxide [20 µg/mL]	Pyridaben [2 µg/mL]
Pyrethrins [10 µg/mL]	Spiroxamine [4 µg/mL]
Spirotetramat [2 µg/mL]	Trifloxystrobin [2 µg/mL]
Thiamethoxam [2 µg/mL]	

